## Wirosil® / Wirosil® plus

### Dublier-Silikon

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### Instructions for use

#### General information

Wirosil® is an addition-cured two-component silicone – Ratio 1 : 1 – that reproduces master models extremely accurately due to its excellent dimensional stability. With economy flask and stabilisation insert it enables work to be carried out easily and reliably without wasting material. Ideal for duplicating milled areas in combination work.

#### Processing

Wirosil®: Cut off the nozzles above the marking.

**After use, always close the nozzles with the caps.**

1. Mix the Wirosil® 1 and 2 components in the measuring and mixing cup in a ratio of 1:1 until showing a homogeneous colour. No colour reams should be seen. Observe the workable time of about 5 minutes (Wirosil®) and approx. 3 min 30 sec (Wirosil®plus).

2. Pour very gently, in a thin stream, the mixed Wirosil® between the master model and the flask border. Follow the flask manufacturer's instructions.

3. Depending on the room-temperature, Wirosil® is hardened after 30 to 40 minutes. It is advantageous to harden the Wirosil®-form at about 4 bar in the pressure unit.

Depending on the room-temperature, Wirosil®plus is hardened after 10 to 12 minutes. It is advantageous to harden the Wirosil®plus-form at about 4 bar in the pressure unit.

**Important:** When working with a pressure unit it has to be observed that the Wirosil®-form and the duplicating model are made under the same pressure-conditions.

4. Blow with compressed air (blow-off nozzle) between the **Wirosil®** mould and the master model. **Then draw out the master model.**

**Note with stronger deformations:** For obtaining the full resilience of duplicating silicone, the duplicating should be filled with investmentmaterial about 30 min after removing the master model.

5. For investment material spray Aurofilm wetting agent onto the Wirosil® model. Let the excess Aurofilm liquid dry at a weak air pressure (blow-off nozzle).

**Important:** With this technique Aurofilm has to be dried completely, as otherwise the investment material would not wet sufficiently the Wirosil® mould: formation of bubbles!

6. Mix the cobalt-chrome investment material, plaster or resin for models in accordance with processing instructions and pour it slowly into the Wirosil® mould standing on the vibrator.

**Important:** When a pressure pot is used, the same pressure conditions must be applied when making the Wirosil® model and the duplicate model.

7. After setting of the investment material blow with compressed air (blow-off nozzle), between the Wirosil® model and the duplicate model. Then draw out the duplicate model.

#### Hints for saving material

1. Grind master model base to be at least 1 cm high!

2. Put small Wirosil® rests remaining from previous duplications in the base such as to leave sufficient space for the liquid Wirosil®.

3. Pour only a quantity of Wirosil® into the flask sufficient to cover by 2 mm the upper edge of the teeth row of the master model.

For duplicating the **Wirosil® duplicating flask system** is highly recommended which enables effective and economical work.

#### Availability

#### Basic assortment **REF** 52000

1 bottle (1 kg) Wirosil® 1

1 bottle (1 kg) Wirosil® 2

1 Wirosil® measuring and mixing bowl

1 Wirosil® spatula

1 Wirosil® duplicating flask, large

1 Wirosil® duplicating flask, small

1 atomizing bottle (100 ml) wetting agent Aurofilm

1 atomizing bottle (100 ml) model spray Durofluid

1 Wirosil® processing instructions

1 instruction for use of duplicating flasks

#### Individual package

Wirosil® duplicating silicone – **REF** 52001

1 bottle (1 kg) Wirosil® 1 white

1 bottle (1 kg) Wirosil® 2 light blue

Wirosil®plus **REF** 54854

1 bottle (1 kg) Wirosil®plus 1 white

1 bottle (1 kg) Wirosil®plus 2 light blue

#### Large package

Wirosil® duplicating silicone – **REF** 51995

1 can (10 kg) Wirosil® 1

1 can (10 kg) Wirosil® 2

#### Accessories

1 Wirosil® duplicating flask, small **REF** 52072

1 Wirosil® duplicating flask, large **REF** 52083

1 stabilizing ring with 3 inserts for small duplicating flask **REF** 52079

large duplicating flask **REF** 52084

1 atomizing bottle (100 ml) Aurofilm wetting agent **REF** 52019

1 atomizing bottle (100 ml) Durofluid model spray **REF** 52008

#### Wirosil®-data

Processing time	about 5 min
Mixing time	about 30 sec
Setting time (22 °C or 71 °F)	30–40 min
Shore-A hardness (1 hour)	17
Ductile yield	250 <span> </span> %
Contraction (DIN 14356)	0.03 <span> </span> %

#### Wirosil®plus-data

Processing time	about 3 min 30 sec
Mixing time	about 30 sec
Setting time (22 °C)	10 – 12 min
Shore-A-hardness (1 hour)	20
Recovery from deformation	99,8 <span> </span> %
Contraction (DIN 14356)	0,01 <span> </span> %

#### Warranty

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values.

Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.

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