

**1 Marker (Dreikantbohrer)**

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

**1 Marker (initial drill)**


Mark the bone and carry out initial drilling for directional orientation. Check with parallel post.

**1a Gap cutter disc\* (optional with spreading)**


Gap cutter disc\* for initial splitting of narrow alveolar ridges.

**1b Gap cutter burr\* (optional with spreading)**

Gap cutter burr\* for deep splitting of alveolar ridges.


**2 Final depth drill**

Final enlargement of the cavity (when using first-generation Mini-/Osseo<sup>Plus</sup> Trays the shape drill following in the tray must also be used!).


\*Gap cutter tools are not included in the current Mini-/Osseo<sup>Plus</sup> Tray and must be purchased separately!

**3 Threadformer (to condense or spread the bone\*\*)**


Form (condense) a thread down to the floor of the cavity or spread a narrow ridge (threadformers with successively increasing diameters can be used until the desired primary stability or distension is achieved\*\*). Can be performed manually or with a handpiece.

**4 Implant insertion – manually or with a handpiece**

Insert the implant into the bone with the ratchet or insertion tool used with a handpiece (machined section remains above the bone). Torque ≤ 40 Ncm.

**5 Placement of the cover screw**

Pick up the cover screw with the hexagon screwdriver, insert into the implant and hand-tighten.

\*\*If the desired primary stability/distension is not achieved after using the Mini-Line threadformers, the threadformers of the BEGO Semados® RI-Line (not included in the current Mini-/Osseo<sup>Plus</sup> tray!) can be optionally used in ascending order until the desired result is achieved in order, for example, to place a BEGO Semados® RI implant with a larger diameter!

**BEGO Semados® Mini-Line**

Indication	Diameter
• Single tooth replacement in the maxilla (lateral incisor)	— Ø 2.9 / 3.1 mm
• Single tooth replacement in the mandible (incisors)	— Ø 2.9 / 3.1 mm
• Bar restoration maxilla / mandible	— Ø 2.7 / 2.9 / 3.1 mm (terminal implants Ø 3.1 mm)
• Ball attachment restoration maxilla / mandible	— Ø 2.9 / 3.1 mm
• Easy-Con restoration maxilla / mandible	— Ø 3.1 mm

**Drill sequence**


Implant	Marker	Gap cutter (optional with spreading)	Depth drill	Shape drill (only in first-generation Mini-/Osseo <sup>Plus</sup> trays!)	Threadformer
Mini 2.7	x	(x)	x	x	x
Mini 2.9	x	(x)	x	x	x
Mini 3.1	x	(x)	x	x	x

x = required (x) = optional

**Speeds**

Article	Speed [rpm]
• Depth drill / Shape drill	— 800
• Threadformer (non-ablative)	— 15-20
• Implant insertion	— 15

**BEGO Implant Systems GmbH & Co. KG**  
 Wilhelm-Herbst-Str. 1  
 28359 Bremen, Germany  
 Tel. +49 421 2028-246  
 Fax +49 421 2028-265  
 E-mail info@bego-implantology.com  
 www.bego.com

**The preferred system for narrow ridges and gaps!**

- Stable implant-abutment connection**  
Conical, external hex connection enables correction of implant angulations
- Microthread in the neck region**  
For improved load transfer in the crestal bone
- Thread (bionic design)**  
Lateral bone compression to increase primary stability
- Surface**  
High-purity and homogenous TiPure<sup>Plus</sup> surface

## THE BEGO SEMADOS® MINI IMPLANT

### An effective implant system for narrow ridge restorations and restorations in narrow spaces

BEGO Semados® Mini implants: The last word in modern and cost-effective dental implantology!

- Cost-effective
- Simple to use
- Practical
- Quality product manufactured entirely in Germany
- Clinically tested
- Particularly recommended for narrow ridge restorations and restorations in narrow spaces
- Bone-friendly bionic thread design

#### The BEGO Semados® Mini implant at a glance

- Made from Grade 4 commercially pure titanium
- Root-shaped
- Diameters: 2.7 · 2.9 · 3.1 mm
- Lengths: 11.5 · 13 · 15 mm
- Excellent osseointegration
- High-purity, homogenous TiPure<sup>Plus</sup> surface
- Bone spreading and bone condensing to enhance primary stability
- Bionic design: Stress reduction for the bone and increase in implant stability

[www.bego.com](http://www.bego.com)

#### The special implant system

- Shorter treatment time thanks to avoidance of time-consuming augmentation
- Low-cost form of treatment
- Permanent restoration of edentulous, atrophied alveolar ridges with prefabricated bar system
- Single tooth and ball attachment restorations possible with implant diameters of 2.9 and 3.1 mm
- Easy-Con Mini restoration with implant diameter of 3.1 mm
- Large selection of prosthetic solutions available
- Stress-free CAD/CAM bar constructions available

Partners in Progress

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REF 84337/02 - CP - © 2017 by BEGO - 2017-01