

Interface Overview*



CADAbut F



CADAbut D



CADBase
implant
niveau



CADBase
abutment
niveau

Implant diameter	Prosthetic platform	BEGO Titan Grade 5	Wirobond® MI+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
BEGO Implant Systems – Semados® S/RI-Line – Sub-Tec (without Platform Switch)														
3,25; 3,75	3,67													
4,1	4,1	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,5	4,5													
5,5	5,5													
BEGO Implant Systems – Semados® SC/SCX/RS/R SX/RI – PS-Line (with Platform Switch)														
3,0	3,0													
3,25; 3,75	3,67													
4,1	4,1	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,5	4,5													
5,5	5,5													
BEGO Implant Systems – Semados® Mini (only for bars)														
2,7; 2,9; 3,1	3,5											√ ^{1) 2)}	√ ^{1) 2)}	
BEGO Implant Systems – Semados® MultiPlus														
3,25; 3,75; 4,1; 4,5; 5,5	abutment niveau												√ ^{1) 2)}	√ ^{1) 2)}
Biocomp (Range 10.1) Implant Level/Abutment Level														
3,4; 4,0	4,0											√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,6; 5,2	5,2													√ ^{1) 2)}
BIOMET 3i (Range 4.1) Osseotite Certain/NanoTite Certain														
3,25; 3/4/3; 4/3	3,4													
4,0; 4/5/4; 5/4	4,1					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
5,0; 5/6/5; 6/5	5,0	√ ^{1) 2)}												
6,0; 5/6	6,0													
BIOMET 3i (Range 4.2) Osseotite/NanoTite														
3,25	3,4													
3,75; 4,0; 3/4	4,1	√ ^{1) 2)}				√ ¹⁾	√ ¹⁾	√ ¹⁾	√ ¹⁾		√ ^{1) 2)}	√ ^{1) 2)}		
5,0; 4/5	5,0													
CAMLOG (Range 8.1) CAMLOG SCREW-LINE/ROOT-LINE														
3,3	3,3													
3,8	3,8													
4,3	4,3	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
5,0	5,0													
6,0	6,0													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.

¹⁾ 3Shape ²⁾ ecoxad

BEGO INTERFACE OVERVIEW

Interface Overview*



CADAbut F



CADAbut D



CADBase
implant
niveau

CADBase
abutment
niveau

Implant diameter	Prosthetic platform	BEGO Titan Grade 5	Wirobond® MI+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
CAMLOG (Range 8.2) CONELOG SCREW-LINE														
3,3	3,3													
3,8	3,8					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}					
4,3	4,3													
5,0	5,0													
CAMLOG (Range 8.3) iSy														
3,8; 4,4; 5,0						√ ¹⁾	√ ¹⁾	√ ¹⁾						
DENTSPLY Implants (Range 2.1) FRIALIT/XiVE														
3,4	3,4													
3,8	3,8					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
4,5	4,5													
5,5	5,5													
DENTSPLY Implants (Range 2.2) Ankylos														
3,5; 4,5; 5,5; 7,0							√ ¹⁾	√ ¹⁾	√ ¹⁾				√ ^{1) 2)}	√ ^{1) 2)}
DENTSPLY Implants (Range 5.1) ASTRA TECH OsseoSpeed														
3,5; 4,0	3,5/4,0					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
4,5; 5,0	4,5/5,0													
DENTSPLY Implants (Range 5.2) ASTRA TECH OsseoSpeed EV														
3,0	3,0													
3,6	3,6													
4,2	4,2					√ ^{1) 2)}							√ ^{1) 2)}	√ ^{1) 2)}
4,8	4,8													
5,4	5,4													
Dyna Dental (Range 9.1) Dyna Helix														
3,2; 3,6; 4,2; 5,0	3,6/4,0										√ ^{1) 2)}	√ ^{1) 2)}		
Neoss (Range 11.1) ProActive														
3,5; 4,0; 4,5; 5,0; 5,5	4,0					√ ^{1) 2)}					√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.1) NobelActive														
3,5	NP					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,3; 5,0	RP													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecoxad

BEGO INTERFACE OVERVIEW

Interface Overview*



CADAbut F



CADAbut D



CADBase
implant
niveau

CADBase
abutment
niveau

Implant diameter	Prosthetic platform	BEGO Titan Grade 5	Wirobond® MI+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
Nobel Biocare (Range 3.2) NobelReplace/Replace Select/NobelSpeedy Replace														
3,5	NP													
4,0; 4,3	RP	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
5,0	WP													
6,0	6,0													
Nobel Biocare (Range 3.3) Brånemark System/NobelSpeedy Shorty/Groovy														
3,3	NP													
3,75; 4,0	RP	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
5,0	WP													
Nobel Biocare (Range 3.4) Multi-unit Abutment														
Nobel Active														
Nobel Replace														
Brånemark System	abutment niveau												√ ^{1) 2)}	√ ^{1) 2)}
ASTRA TECH OsseoSpeed														
Straumann synOcta														
Straumann BoneLevel														
Phibo (Range 7.1) TSH														
S2 3,3	3,3													
S3 3,6; S4 4,2	4,0	√ ^{1) 2)}									√ ^{1) 2)}	√ ^{1) 2)}		
S5 4,8	5,0													
Straumann (Range 1.1) synOcta/Standard/Standard Plus/Tapered Effect														
3,3; 4,1; 4,8	RN 4,8	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,8	WN 6,5													
Straumann (Range 1.2) Straumann BoneLevel														
3,3	NC	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
4,1; 4,8	RC													
Straumann (Range 1.3) Straumann Narrow Neck CrossFit														
3,5	NNC	√ ^{1) 2)}				√ ¹⁾	√ ¹⁾	√ ¹⁾						
Zimmer Dental (Range 6.1) Tapered Screw-Vent														
3,7; 4,1	3,5													
4,7	4,5	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
6,0	5,7													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecoxad