

Interface Overview*



CADAbut F



CADAbut D



CADBase
implant
niveau







CADBase
abutment
niveau

Implant diameter	BEGO Titan Grade 5	Wirobond® MI+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
BEGO Implant Systems – Semados® SC/SCX/RS/R SX/RI – PS-Line (with Platform Switch)													
3,0													
3,25; 3,75													
4,1	√ ^{1) 2)}												
4,5			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
5,5													
BEGO Implant Systems – Semados® S/RI-Line – Sub-Tec (without Platform Switch)													
3,25; 3,75													
4,1	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}										
4,5	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
5,5													
BEGO Implant Systems – Semados® Mini (only for bars)													
2,7; 2,9; 3,1												√ ^{1) 2)}	√ ^{1) 2)}
BEGO Implant Systems – Semados® MultiPlus													
3,25; 3,75; 4,1; 4,5; 5,5 abutment niveau													√ ^{1) 2)} √ ^{1) 2)}
Biocomp (Range 10.1) Implant Level/Abutment Level													
3,4; 4,0												√ ^{1) 2)}	√ ^{1) 2)}
4,6; 5,2												√ ^{1) 2)}	√ ^{1) 2)}
BIOMET 3i (Range 4.1) Osseotite Certain/NanoTite Certain													
3,25; 3/4/3; 4/3													
4,0; 4/5/4; 5/4					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}
5,0; 5/6/5; 6/5	√ ^{1) 2)}												
6,0; 5/6													
BIOMET 3i (Range 4.2) Osseotite/NanoTite													
3,25													
3,75; 4,0; 3/4	√ ^{1) 2)}				√ ¹⁾	√ ¹⁾	√ ¹⁾	√ ¹⁾				√ ^{1) 2)}	√ ^{1) 2)}
5,0; 4/5													
Bredent medical (Range 15.1) Sky													
3,5													
4,0; 4,5; 5,5	√ ^{1) 2)}												

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecocad

BEGO INTERFACE OVERVIEW

Interface Overview*													
	 CADAbut F			 CADAbut D				 CADBase implant niveau		 CADBase abutment niveau			
Implant diameter	BEGO Titan Grade 5	Wirobond® Ml+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
CAMLOG (Range 8.1) CAMLOG SCREW-LINE /ROOT-LINE													
3,3													
3,8													
4,3	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
5,0													
6,0													
CAMLOG (Range 8.2) CONELOG SCREW-LINE													
3,3													
3,8													
4,3					√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}					
5,0													
CAMLOG (Range 8.3) iSy													
3,8; 4,4; 5,0					√ ¹⁾	√ ¹⁾	√ ¹⁾						
Dentaurum (Range 14.1) tioLogic®													
3,3													
3,7; 4,2					√ ^{1) 2)}								
4,8; 5,5													
DENTSPLY Implants (Range 2.1) FRIALIT/XiVE													
3,4													
3,8	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
4,5													
5,5													
DENTSPLY Implants (Range 2.2) Ankylos													
3,5; 4,5; 5,5; 7,0					√ ¹⁾	√ ¹⁾	√ ¹⁾					√ ^{1) 2)}	√ ^{1) 2)}
DENTSPLY Implants (Range 5.1) ASTRA TECH OsseoSpeed													
3,5; 4,0	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
4,5; 5,0		√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
DENTSPLY Implants (Range 5.2) ASTRA TECH OsseoSpeed EV													
3,0													
3,6													
4,2	√ ^{1) 2)}											√ ^{1) 2)}	√ ^{1) 2)}
4,8													
5,4													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecocad

BEGO INTERFACE OVERVIEW

Interface Overview*



CADAbut F



CADAbut D



CADBase
implant
niveau

CADBase
abutment
niveau

Implant diameter	BEGO Titan Grade 5	Wirobond® MI+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
Dyna Dental (Range 9.1) Dyna Helix													
3,2; 3,6; 4,2; 5,0										√ ^{1) 2)}	√ ^{1) 2)}		
Neoss (Range 11.1) ProActive													
3,5; 4,0; 4,5; 5,0; 5,5	√ ^{1) 2)}									√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.1) NobelActive													
3,5 4,3; 5,0	√ ^{1) 2)}		√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.2) NobelReplace/Replace Select/NobelSpeedy Replace													
3,5 4,0; 4,3 5,0 6,0	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.3) Brånemark System/NobelSpeedy Shorty/Groovy													
3,3 3,75; 4,0 5,0	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.4) Multi-unit Abutment													
Nobel Active, abutment niveau Nobel Replace, abutment niveau Brånemark System, abutment niveau ASTRA TECH OsseoSpeed, abutment niveau Straumann synOcta, abutment niveau Straumann BoneLevel, abutment niveau													√ ^{1) 2)} √ ^{1) 2)}
Phibo (Range 7.1) TSH													
S2 3,3 S3 3,6; S4 4,2 S5 4,8	√ ^{1) 2)}									√ ^{1) 2)}	√ ^{1) 2)}		
Straumann (Range 1.1) synOcta/Standard/Standard Plus/Tapered Effect													
3,3; 4,1; 4,8 4,8	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecoxad

BEGO INTERFACE OVERVIEW

Interface Overview*



CADAbut F



CADAbut D



CADBase
implant
niveau

CADBase
abutment
niveau

Implant diameter	BEGO Titan Grade 5	Wirobond® MI+	BeCe® CAD Zirkon XH	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
Straumann (Range 1.2) Straumann BoneLevel													
3,3 4,1; 4,8	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
Straumann (Range 1.3) Straumann Narrow Neck CrossFit													
3,5	√ ^{1) 2)}				√ ¹⁾	√ ¹⁾	√ ¹⁾						
Zimmer Dental (Range 6.1) Tapered Screw-Vent													
3,7; 4,1 4,7 6,0	√ ^{1) 2)}				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecoxad

BEGO Medical GmbH

Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany · Hotline +49 421 2028-200

Always want to have a look at all of BEGO's news?

Sign up for our Newsletter here: www.bego.com/newsletter

