

Interface Overview*													
		CADAbut Full		CADAbut Duo				CADBase implant niveau		CADBase abutment niveau			
Implant diameter	Color coding	BEGO Titan Grade 5	Wirobond® M+	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
BEGO Implant Systems – Semados® SC/SCX/RS/R SX/RI – PS-Line (with Platform Switch)													
RS/R SX 3,0	■												
SC/SCX 3,25	■												
SC/SCX/RS/R SX/RI 3,25; 3,75	■	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
SC/SCX/RS/R SX/RI 4,1	■	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
SC/SCX/RS/R SX/RI 4,5	■												
SC/SCX/RS/R SX/RI 5,5	■												
BEGO Implant Systems – Semados® S/RI-Line – Sub-Tec (without Platform Switch)													
S 3,25	■												
S/RI 3,75	■												
S/RI 4,1	■	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
S/RI 4,5	■												
S/RI 5,5	■												
BEGO Implant Systems – Semados® Mini (only for bars)													
2,7; 2,9; 3,1	■									√ ^{1) 2)}	√ ^{1) 2)}		
BEGO Implant Systems – Semados® Multi^{Plus}													
3,25; 3,75; 4,1; 4,5; 5,5 Abutment Niveau												√ ^{1) 2)}	√ ^{1) 2)}
Biocomp (Range 10.1) Implant Level/Abutment Level													
3,4; 4,0										√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,6; 5,2													
BIOMET 3i (Range 4.1) Osseotite Certain/NanoTite Certain													
3,25; 3/4/3; 4/3	■												
4,0; 4/5/4; 5/4	■	√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
5,0; 5/6/5; 6/5	■												
6,0; 5/6	■												
BIOMET 3i (Range 4.2) Osseotite/NanoTite													
3,25	■												
3,75; 4,0; 3/4	■	√ ^{1) 2)}			√ ¹⁾	√ ¹⁾	√ ¹⁾	√ ¹⁾		√ ^{1) 2)}	√ ^{1) 2)}		
5,0; 4/5	■												
Bredent medical (Range 15.1) Sky													
3,5		√ ^{1) 2)}											
4,0; 4,5; 5,5													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.
¹⁾ 3Shape ²⁾ ecocad

BEGO INTERFACE OVERVIEW

Interface Overview*													
		CADAbut Full		CADAbut Duo				CADBase implant niveau		CADBase abutment niveau			
Implant diameter	Color coding	BEGO Titan Grade 5	Wirobond® M1+	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
CAMLOG (Range 8.1) CAMLOG SCREW-LINE/ROOT-LINE													
3,3	■												
3,8	■												
4,3	■	√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
5,0	■												
6,0	■												
CAMLOG (Range 8.2) CONELOG SCREW-LINE													
3,3	■												
3,8	■												
4,3	■				√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}					
5,0	■												
CAMLOG (Range 8.3) iSy													
3,8; 4,4; 5,0					√ ¹⁾	√ ¹⁾	√ ¹⁾						
Dentaurum (Range 14.1) tioLogic													
3,3													
3,7; 4,2		√ ^{1) 2)}											
4,8; 5,5													
DENTSPLY Implants (Range 2.1) FRIALIT/XiVE													
3,4	■												
3,8	■												
4,5	■	√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
5,5	■												
DENTSPLY Implants (Range 2.2) Ankylos													
3,5; 4,5; 5,5; 7,0					√ ¹⁾	√ ¹⁾	√ ¹⁾					√ ^{1) 2)}	√ ^{1) 2)}
DENTSPLY Implants (Range 5.1) ASTRA TECH OsseoSpeed													
3,5; 4,0	■												
4,5; 5,0	■	√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.

¹⁾ 3Shape ²⁾ ecocad

BEGO INTERFACE OVERVIEW

Interface Overview*													
		CADAbut Full		CADAbut Duo				CADBase implant niveau		CADBase abutment niveau			
Implant diameter	Color coding	BEGO Titan Grade 5	Wirobond® M+	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
DENTSPLY Implants (Range 5.2) ASTRA TECH OsseoSpeed EV													
3,0													
3,6													
4,2				√ ^{1) 2)}								√ ^{1) 2)}	√ ^{1) 2)}
4,8													
5,4													
Dyna Dental (Range 9.1) Dyna Helix													
3,2; 3,6; 4,2; 5,0										√ ^{1) 2)}	√ ^{1) 2)}		
Neoss (Range 11.1) ProActive													
3,5; 4,0; 4,5; 5,0; 5,5			√ ^{1) 2)}							√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.1) NobelActive													
3,5													
4,3; 5,0			√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
Nobel Biocare (Range 3.2) NobelReplace/Replace Select/NobelSpeedy Replace													
3,5													
4,0; 4,3													
5,0			√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
6,0													
Nobel Biocare (Range 3.3) Brånemark System/NobelSpeedy Shorty/Groovy													
3,3													
3,75; 4,0			√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
5,0													
Nobel Biocare (Range 3.4) Multi-unit Abutment													
Nobel Active, Abutment Niveau													
Nobel Replace, Abutment Niveau													
Brånemark System, Abutment Niveau													
ASTRA TECH OsseoSpeed, Abutment Niveau												√ ^{1) 2)}	√ ^{1) 2)}
Straumann synOcta, Abutment Niveau													
Straumann BoneLevel, Abutment Niveau													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.

¹⁾ 3Shape ²⁾ ecocad

BEGO INTERFACE OVERVIEW

Interface Overview*													
		CADAbut Full		CADAbut Duo				CADBase implant niveau		CADBase abutment niveau			
Implant diameter	Color coding	BEGO Titan Grade 5	Wirobond® Ml+	IPS e.max CAD	BEGO Zirkon LT	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT bridge bonded to titanium no Hex base	BEGO Titan Grade 5	Wirobond® M+	BEGO Titan Grade 5	Wirobond® M+
Phibo (Range 7.1) TSH													
S2 3,3													
S3 3,6; S4 4,2		√ ^{1) 2)}								√ ^{1) 2)}	√ ^{1) 2)}		
S5 4,8													
Straumann (Range 1.1) synOcta/Standard/Standard Plus/Tapered Effect													
3,3; 4,1; 4,8		√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}
4,8													
Straumann (Range 1.2) Straumann BoneLevel													
3,3		√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}		√ ^{1) 2)}	√ ^{1) 2)}		
4,1; 4,8													
Straumann (Range 1.3) Straumann Narrow Neck CrossFit													
3,5		√ ^{1) 2)}			√ ¹⁾	√ ¹⁾	√ ¹⁾						
Zimmer Dental (Range 6.1) Tapered Screw-Vent													
3,7; 4,1										√ ^{1) 2)}	√ ^{1) 2)}		
4,7		√ ^{1) 2)}			√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}	√ ^{1) 2)}					
6,0													

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems.

¹⁾ 3Shape ²⁾ ecoxad

BEGO Medical GmbH

Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany
 Hotline +49 421 2028-200
 Tel. +49 421 2028-0 · Fax +49 421 2028-174
 Email info@bego-medical.com · www.bego.com

Always want to have a look at all of BEGO's news?
 Sign up for our Newsletter here: www.bego.com/newsletter

