

Parallelisierungsposten für Mini-Line
Art.-Nr. 56972

Dreikantbohrer
Art.-Nr. 55714

DD (Tiefenbohrer) Mini 2,7
Art.-Nr. 57995

Gewindeformer Mini 2,7
Art.-Nr. 56979

Eindrehwerkzeug Implantate, kurz, mit Ratschenanschluss
Art.-Nr. 55831

Eindrehwerkzeug Implantate, lang, mit Ratschenanschluss
Art.-Nr. 55830

Eindrehwerkzeug Implantate, lang, mit Winkelstückanschluss
Art.-Nr. 55832

Bohrerverlängerung einteilig, ohne Innenkühlung
Art.-Nr. 57265

DD (Tiefenbohrer) Mini 2,9
Art.-Nr. 57996

Gewindeformer Mini 2,9
Art.-Nr. 56980

Sechskantschlüssel 1,25 mm, kurz, mit Ratschenanschluss
Art.-Nr. 53949

Sechskantschlüssel 1,25 mm, lang, mit Ratschenanschluss
Art.-Nr. 53948

Sechskantschlüssel 1,25 mm, kurz, mit Winkelstückanschluss
Art.-Nr. 55563

Sechskantschlüssel 1,25 mm, lang, mit Winkelstückanschluss
Art.-Nr. 55564

Drilling Depth
11.5, 15, 13

CE

BEGO Semados® Mini-Line

SEMADOS

OUT

Ratsche/Drehmomentschlüssel
Art.-Nr. 55799

Führungsschlüssel*
Art.-Nr. 53897

DD (Tiefenbohrer) Mini 3,1
Art.-Nr. 57997

Gewindeformer Mini 3,1
Art.-Nr. 56981

* Instrument im Container, unterhalb des sichtbaren Einsatzes.

Paralleling post for Mini-Line
Art. no. 56972

Initial drill
Art. no. 55714

DD (Depth drill) Mini 2.7
Art. no. 57995

Threadformer Mini 2.7
Art. no. 56979

Insertion tool for implants, short, with ratchet connection
Art. no. 55831

Insertion tool for implants, long, with ratchet connection
Art. no. 55830

Insertion tool for implants, long, with contra-angle connection
Art. no. 55832

Drill extension, one-piece, without internal cooling
Art. no. 57265

Hexagon screw-driver 1.25 mm, short, with ratchet connection
Art. no. 53949

Hexagon screw-driver 1.25 mm, long, with ratchet connection
Art. no. 53948

Hexagon screw-driver 1.25 mm, short, with contra-angle connection
Art. no. 55563

Hexagon screw-driver 1.25 mm, long, with contra-angle connection
Art. no. 55564

DD (Depth drill) Mini 2.9
Art. no. 57996

Threadformer Mini 2.9
Art. no. 56980

DD (Depth drill) Mini 3.1
Art. no. 57997

Threadformer Mini 3.1
Art. no. 56981

Ratchet/Torque spanner
Art. no. 55799

Guiding wrench*
Art. no. 53897

Drilling Depth
11.5
15
13

CE

BEGO Semados® Mini-Line

SEMADOS

OUT

Hexagon Screwdriver

Depth Drill

Thread-former

2.7

2.9

3.1

Pilot Marker

Parallel Pin

Insertion Tool

Drill Extension

Add-on

* Instrument in the container, below the visible tray.