


The evolution of the BEGO Semados® S-Line



BEGO SEMADOS® S IMPLANTS HAVE BECOME SC AND SCX IMPLANTS.

Your trusted system with an improved design

BEGO Semados® SC/SCX implants: The last word in modern and cost-effective dental implantology

- Quick & functional
- Reliable & trusted
- Simple & intuitive
- Universal & successful
- Varied & sophisticated
- 100% German quality product

The BEGO Semados® SC/SCX implant at a glance

- Made of pure grade 4 titanium
- Diameter: 3.25 · 3.75 · 4.1 · 4.5 · 5.5 mm
- Lengths: 7 · 8.5 · 10 · 11.5 · 13 · 15 mm
- Diameter 5.5 in length 7 mm to prevent vertical augmentations
- High purity, homogeneous TiPure^{Plus} surface
- Cylindrical implant with rounded apex to protect anatomical structures

An overview of the SC/SCX implant prosthetics connection

- Platform Switch
- Established conical internal connection with 45° conical taper and hexagon socket for anti-rotation protection
- Inner cone ensures tight fit of the restoration
- CAD/CAM customised prosthetic components made of various materials

www.bego.com/scx

Partners in Progress



Addresses · Delivery address · Invoice address · Customer number (if available)

I hereby order the products listed below.
Please send the delivery to the above address.

BEGO Semados® SC Implants			
Art. No.	Colour code	Product	Quantity
BEGO Semados® SC implants 3.25			
58160	■	SC 3.25 L 8.5	
58161	■	SC 3.25 L 10	
58162	■	SC 3.25 L 11.5	
58163	■	SC 3.25 L 13	
58164	■	SC 3.25 L 15	
BEGO Semados® SC implants 3.75			
58165	■	SC 3.75 L 8.5	
58166	■	SC 3.75 L 10	
58167	■	SC 3.75 L 11.5	
58168	■	SC 3.75 L 13	
58169	■	SC 3.75 L 15	
BEGO Semados® SC implants 4.1			
58170	■	SC 4.1 L 7	
58171	■	SC 4.1 L 8.5	
58172	■	SC 4.1 L 10	
58173	■	SC 4.1 L 11.5	
58174	■	SC 4.1 L 13	
58175	■	SC 4.1 L 15	
BEGO Semados® SC implants 4.5			
58176	■	SC 4.5 L 7	
58177	■	SC 4.5 L 8.5	
58178	■	SC 4.5 L 10	
58179	■	SC 4.5 L 11.5	
58180	■	SC 4.5 L 13	
58181	■	SC 4.5 L 15	
BEGO Semados® SC implants 5.5			
58182	■	SC 5.5 L 7	
58183	■	SC 5.5 L 8.5	
58184	■	SC 5.5 L 10	
58185	■	SC 5.5 L 11.5	
58186	■	SC 5.5 L 13	

BEGO Semados® SCX implants			
Art. No.	Colour code	Product	Quantity
BEGO Semados® SCX implants 3.25			
58260	■	SCX 3.25 L 8.5	
58261	■	SCX 3.25 L 10	
58262	■	SCX 3.25 L 11.5	
58263	■	SCX 3.25 L 13	
58264	■	SCX 3.25 L 15	
BEGO Semados® SCX implants 3.75			
58265	■	SCX 3.75 L 8.5	
58266	■	SCX 3.75 L 10	
58267	■	SCX 3.75 L 11.5	
58268	■	SCX 3.75 L 13	
58269	■	SCX 3.75 L 15	
BEGO Semados® SCX implants 4.1			
58270	■	SCX 4.1 L 7	
58271	■	SCX 4.1 L 8.5	
58272	■	SCX 4.1 L 10	
58273	■	SCX 4.1 L 11.5	
58274	■	SCX 4.1 L 13	
58275	■	SCX 4.1 L 15	
BEGO Semados® SCX implants 4.5			
58276	■	SCX 4.5 L 7	
58277	■	SCX 4.5 L 8.5	
58278	■	SCX 4.5 L 10	
58279	■	SCX 4.5 L 11.5	
58280	■	SCX 4.5 L 13	
58281	■	SCX 4.5 L 15	
BEGO Semados® SCX implants 5.5			
58282	■	SCX 5.5 L 7	
58283	■	SCX 5.5 L 8.5	
58284	■	SCX 5.5 L 10	
58285	■	SCX 5.5 L 11.5	
58286	■	SCX 5.5 L 13	

Place, date, signature

BEGO Implant Systems GmbH & Co. KG
 Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany
 Tel. +49 421 2028-246 · Fax +49 421 2028-265
 E-Mail info@bego-implantology.com · www.bego.com