



CAD / CAM Solutions

# BEGO SCANNERS AND SOFTWARE

Valid from  
October 2026

Partners in Progress



Dear Dental Experts,

Digital solutions are an integral part of dentistry and dental technology today. What was once the technology of tomorrow is state-of-the-art today. We no longer talk about whether a new digital system is as good as the classic analogue technology, but rather what additional benefits and possibilities it offers.

Next-generation laboratory scanners are even more focussed than before on increasing efficiency. With CAD software, the number of possible indications is constantly increasing. Ways to achieve your final design within the software are increasing at the same time as efficient shorter design flows are being offered. With intraoral scans, we no longer talk about whether they are as accurate as conventional impressions, but rather what new value creation opportunities arise from this. And last but not least, we are talking about the precise preoperative planning of implant positions in a ‘forward thinking’ workflow in which we look ahead together – towards the perfect prosthetic restoration.

Digital workflows are therefore the key to even closer integration between dentists and dental technicians. When data records and planning proposals can be exchanged within seconds, treatments can be discussed independently of location and you take patient care to a new level.

At BEGO, we are continuously working on further developing these workflows and providing you with the best possible digital support. For this, we offer you an ideal basis with our digital portfolio of scanners and software. Convince yourself!

Best regards

**Torsten Bahr**  
Senior Product Manager Digital Platforms

WWW.BEGO.COM

# Contents

## 1

### Medit\* Intraoral Scanners

- 6 The Medit i900 Family
- 14 Medit i700/i700w
- 18 Accessories Medit Scanner

## 2

### Medit Laboratory Scanner and exocad\*

- 22 Laboratory Scanner and CAD Software
- 23 Medit T310, T510 & T710
- 25 exocad Software Packages
- 26 Optional exocad Software Modules

## 3

### 3Shape\* Laboratory Scanner

- 29 Laboratory Scanner and CAD Software
- 30 Choose Your 3Shape Red Scanner – E2, E3, E4, F8
- 31 Choose Your 3Shape Software Packages
- 31 Choose PC
- 32 Choose Your 3Shape LabCare Service
- 33 Choose Optional 3Shape Software Modules
- 34 Choose Optional 3Shape Accessories

## 4

### Implantat and Chairside Planning

- 38 exocad exoplan
- 39 3Shape Implant Studio
- 40 BEGO Guide Planning Center
- 41 exocad chairsideCAD

## 5

### Service

- 44 BEGOconnectApp
- 44 BEGO Education

\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.

# 1

## Medit\* Intraoral Scanners

\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.

# TRUE MOBILITY FOR EVERYDAY PRACTICE

With the Medit i900M, experience the freedom to scan anytime, anywhere, without compromise

# MEDIT i900M

## TRUE MOBILITY. ANYWHERE.

A new scanning  
experience.

### Scan Button

Instantly start or stop  
scans with responsive  
tactile feedback

### Control Button

Navigate views  
hands-free - pan,  
zoom, and rotate directly  
from the handpiece

### Built-in Wi-Fi 6E

Ultra-fast, low-latency scanning.  
Wire-free, hassle-free

### LED Status Indicators

Status colors for scanning,  
standby, and morer

### Wider FOV

Field of View: 18 x 15 mm

### Reusable Tips

Three sizes: S, M, L



### Robust and Durable Construction

Precision-sealed body design  
minimizes gaps, offering  
water-resistant durability

### Battery Management System Built-in

Up to 1 hour of continuous  
use or 30+ scans per charge.  
Features thermal protection  
and real-time alerts



# MOBILITY IN YOUR HANDS

Medit\* i900M



Experience true mobility in your daily workflow with the Medit i900M. Advanced third-generation optics, intelligent filtering, and wireless performance deliver high-resolution scans in a sleek, lightweight design. The ergonomic handpiece and intuitive operation reduce physical strain and save valuable time, even in complex cases.

### Sharper. Faster. Stronger.

- **13rd-Generation 10-Bit Imaging**  
Ultra-realistic images with enhanced color depth and texture detail
- **Clear, Accurate Metal Scanning**  
Captures metal surfaces with outstanding clarity & minimal distortion
- **Effortless Deep Scanning**  
Up to 30 mm scan depth enables reliable capture of restorations in hard-to-reach areas
- **Wireless Speed**  
Wi-Fi 6E delivers fast, stable scanning with minimal interference
- **Dual Optical Cameras**  
Provide sharp, high-precision scans with exceptional clarity
- **Versatility**  
Adapts seamlessly to any clinical environment and workflow

**Versatility**  
Use an iPad, PC or the web – whichever suits your hospital environment best. **To access the full range of features, we'd be happy to recommend a suitable laptop.**

### Minimal learning curve – with iPad app (optional)

The i900M's user interface is so intuitive that any team member can scan with confidence – without the need for any training.

- **No complicated manuals**  
Simply open the app and start scanning
- **Instant start-up**  
No booting, no loading screens. Start and scan straight away
- **Optimised UI**  
A clear layout that shows only the most important tools
- **Intuitive for every team member – efficient scanning without specialist knowledge**

### True mobility

The i900M adapts to real-world situations – ANYWHERE


- **Chairside-Scanning**  
With just one device, you can switch seamlessly between workstations
- **Real-time case sharing**  
Transfer scans instantly to other staff members
- **Real-time patient communication**  
Show 3D scans live with Medit Link Express and build trust immediately through transparency and clear visualisations

### Product details Medit\* i900M

|                       |                            |                                                                                                |
|-----------------------|----------------------------|------------------------------------------------------------------------------------------------|
| Scanning Technology   | Capture Rate               | Fastest capture rate in the Medit ecosystem                                                    |
|                       | Imaging Technology         | 3D-in-motion video technology   3D full color streaming capture                                |
|                       | Anti-fogging Technology    | Adaptive anti-fogging                                                                          |
|                       | Camera Grab                | 10 bit                                                                                         |
| Accuracy              | Full-arch                  | 10.9 µm ± 0.98                                                                                 |
| Handpiece             | Dimensions                 | 209 × 42 × 42 mm                                                                               |
|                       | Total Weight               | 319 g (With a medium-sized tip and a battery)                                                  |
| Tip                   | Size of the tips           | L: 36,1 × 34,1 × 90,8 mm; M: 36,1 × 34,1 × 90,4 mm; S: 36,1 × 34,1 × 90,3 mm (sold separately) |
|                       | Field of View (Resolution) | L: 18 × 15 mm (552 × 460 px); M: 14 × 11.5mm (428 × 352 px) S: 10 × 8 mm (308 × 252 px)        |
|                       | Autoclavable               | Up to 150 times<br>121 °C 30 min<br>134 °C 4 min<br>135 °C 10 min                              |
|                       | Material                   | PPSU Polyphenylsulfone (mostly used for baby bottles)                                          |
|                       | 180° Reversible Tip        | No                                                                                             |
|                       | Automatic Tip Detection    | Yes                                                                                            |
| Connectivity          |                            | Wireless                                                                                       |
| Battery Specification |                            | Working hour: 1 h<br>Charging hour: approx. 2 h                                                |
| REF Medit i900M       |                            | 280736                                                                                         |



Detailed  
Product information



\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.

\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.

# ALL IN ONE INNOVATION

The Medit\* i900 features cutting-edge technology for smoother scanning, easier handling and improved scan accuracy

# MEDIT

## i900/i900c

### THE EVOLUTION OF INTRAORAL SCANNING

A smarter, easier  
way to work

#### TouchBand™

Real-time light and haptic feedback and button-free rotation

#### TouchPad

Customizable command actions and mouse-free smart interaction

#### Fast Preheating

Ready to scan instantly

#### Unibody-Design

Housing design without gaps: Experience the unique haptics together with improved hygiene

#### Wand Tip

Multiple tip sizes to customize patient comfort

#### Intelligent Sensor

LED and haptic feedback

#### 3rd-Generation Optical Engine

Advanced 10-bit camera technology

#### Light & Compact

The handpiece weighs only 165 g and is very compact with a length of 233 mm



TouchBand™ and TouchPad



Maximum compactness for convenient scanning



Also available as i900c without TouchBand™ and TouchPad

\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.



# All IN ONE Package for Usability

Medit\* i900/i900c

With a lightweight, ergonomic design, the i900 ensures anyone in your clinic can master scanning with unparalleled ease. An increased field of view captures a wider area with minimal effort, ensuring speed and accuracy in every session. A new unibody design ensures perfect balance and eliminates contamination worries.

## Revolutionary TouchBand™ and TouchPad Functionality\*\*

Level up your experience with a scanner that redefines the scanning process with an intuitive TouchBand™ and TouchPad. Our ground-breaking 360° TouchBand provides easy navigation with haptic feedback. A button-free operation gives you precise control from any angle and allows you to easily start and stop scanning.

## Ergonomic with Smart Interaction

Light. Compact. Portable. Featherlight at 165 grams (5.8 ounces) and built for easy scanning, our compact, portable scanner offers unbridled freedom.

## ALL IN ONE for your comfort

### Enjoy a hassle-free experience right from the start!

The Medit i900 is the complete scanning solution for your practice. Plus, every i900 is backed by direct Medit service and a 3-year warranty, ensuring your peace of mind. The direct service provided by BEGO in combination with Medit ensures that you receive a carefree scanning experience at all times.

## Simple to set up

### Start scanning right out of the box!

#### The Medit i900/i900c:

- Handpiece and holder
- Four wand tips in two sizes
- 2 flexible, hassle-free cables
- 3-year warranty and direct service from Medit and BEGO



## Product details Medit\* i900/i900c

|                     |                            |                                                                                         |
|---------------------|----------------------------|-----------------------------------------------------------------------------------------|
| Scanning Technology | Capture Rate               | Fastest capture rate in the Medit ecosystem                                             |
|                     | Imaging Technology         | 3D-in-motion video technology I 3D full color streaming capture                         |
|                     | Anti-fogging Technology    | Adaptive anti-fogging                                                                   |
|                     | Camera Grab                | 10 bit                                                                                  |
| Accuracy            | Full-arch                  | 10.9 μm ± 0.98                                                                          |
| Handpiece           | Dimensions                 | 223.4 × 36.7 × 19.7 mm                                                                  |
|                     | Total Weight               | 165 g                                                                                   |
| Tip                 | Size                       | L: 26.7 × 17.1 mm; M: 21.8 × 16.3 mm; S: 17.9 × 13.1 mm (sold separately)               |
|                     | Field of View (Resolution) | L: 18 × 15 mm (552 × 460 px); M: 14 × 11.5mm (428 × 352 px) S: 10 × 8 mm (308 × 252 px) |
|                     | Autoclavable               | Up to 150 times<br>121 °C 30 min<br>134 °C 4 min<br>135 °C 10 min                       |
|                     | Material                   | PPSU Polyphenylsulfone (mostly used for baby bottles)                                   |
|                     | 180° Reversible Tip        | No                                                                                      |
|                     | Automatic Tip Detectiong   | Yes                                                                                     |
| Cable               |                            | C Power Delivery: 2 m & 2.5 m                                                           |
| REF Medit i900      |                            | 280733                                                                                  |
| REF Medit i900c     |                            | 280735                                                                                  |
| REF Medit Cart      |                            | 100243                                                                                  |



Detailed  
Product information

Medit i900/i900c  
incl.  
3-year warranty

\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.  
\*\* only i900

\* This symbol is a commercial designation / registered trademark of a company which is not part of the BEGO company group.

# MAGIC MADE EASY

Get the power to take your practice to the next level!

# MEDIT i700 / i700w

PERFORMANCE

Full arch  
accuracy:  
 $10.9 \mu\text{m} \pm 0.98$  –  
up to 70 fps

## Small Tip

(sold separately)  
For easy application in the  
posterior region for patients  
with a small oral cavity

## Wider Scan Area

Enlarged scanning window captures  
large areas at once while maintaining  
the small size of the scanner head

## Reversible Tip

No need to twist your  
wrist for maxillary scans

## Auto Calibration

Automatic tool makes calibration  
a breeze



Auto Calibration Tool  
(sold separately)

## Remote Control

Conveniently check  
scan data from  
the scanner itself  
and prevent cross-  
contamination

## UV-C LED Disinfection

Lower risk of contamination  
inside the handpiece

## Lightweight & Superbly Balanced


For a stable, comfortable use

## Battery

Scan up to 30 cases  
with 1 battery\*

## 60 GHz Wireless Technology

For high-speed wireless connection



Hub


When the scanner's cable disappears, so do restrictions to scan movements. Strengthened impact resistance minimizes damage from accidental drops. Weight balanced design of the i700 wireless provides a stable, comfortable scanning experience.

\* Performance may vary according to usage conditions.



# A perfect Duo – Medit i700/i700w

Wireless or wired – You decide!



### The Best Battery Life Ever

1 hour of continuous scanning, up to 8 hours in standby mode. Battery time can vary depending on the user environment (temperature, signal strength, etc.) and battery condition.



### All Day Scan with Single Battery


Scan all day without worrying about charging. Scan up to 30 cases with 1 battery, up to 90 cases using 3 included batteries.



8 HOURS


### Intelligent power management function

Intelligent power management function automatically switches the device to sleep mode when not in use, standing by up to 8 hours. It takes just seconds to wake the device for scanning.



### Wired power connection

Once you've used up all 3 included batteries, simply connect the power cable and keep using the device.



### Total scan time of 2 mins\*\*

Maxillary + Mandibular + Occlusion

### Product details Medit\* i700/i700w

|                     |                                                                           |                                                                                                                   |
|---------------------|---------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------|
| Scanning Technology | Scanning Frame                                                            | Up to 70 FPS                                                                                                      |
|                     | Imaging Technology                                                        | 3D-in-motion video technology I 3D full color streaming capture                                                   |
|                     | LED Source                                                                | LED                                                                                                               |
|                     | Anti-Fogging Technology                                                   | Adaptive anti-fogging                                                                                             |
| Accuracy            | Full-Arch                                                                 | 10.9 μm ± 0.98                                                                                                    |
| Handpiece           | Dimensions                                                                | i700: 248 × 44 × 47.4 mm<br>i700w: 313 × 44 × 47.4 mm                                                             |
|                     | Total Weight                                                              | i700: 245 g (Includes tip)<br>i700w: 328 g (Includes battery, tip)                                                |
| Tip                 | Tip Size                                                                  | 22.5 × 17.1 mm I 19.0 × 14.8 mm (Small Tip)                                                                       |
|                     | Scan Area                                                                 | 15 × 13 mm                                                                                                        |
|                     | Autoclavable                                                              | Up to 150 times<br>121 °C 30 min (Gravity Type)<br>134 °C 4 min (Pre-Vacuum Type)<br>135 °C 10 min (Gravity Type) |
|                     | 180° Reversible Tip                                                       | Yes                                                                                                               |
| Special Mode        | Remote Control Mode                                                       | Yes                                                                                                               |
|                     | UV-C Disinfection                                                         | Yes                                                                                                               |
| Connectivity        | i700: Wire, Length 2.0 m<br>i700w: Wireless, USB 3.1 Gen1 (Wireless Hub ) |                                                                                                                   |
| REF                 | i700: 280731<br>i700w: 280730                                             |                                                                                                                   |



Detailed product information

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
\*\* Performance may vary according to usage conditions.

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

Accessories Medit i900/i900c/i900M

| Article                                                      | i900 | i900c | i900M | REF    |
|--------------------------------------------------------------|------|-------|-------|--------|
| Medit i900/i900c/i900M Handpiece Cover                       | ✓    | ✓     | ✓     | 100233 |
| Medit i900/i900c/i900M Large and Medium Tip Set (each 2 pcs) | ✓    | ✓     | ✓     | 100234 |
| Medit i900/i900c/i900M Calibration Tool                      | ✓    | ✓     | ✓     | 100235 |
| Medit i900/i900c Desktop Cradle                              | ✓    | ✓     |       | 100236 |
| Medit i900/i900c Wall Mount Holder                           | ✓    | ✓     |       | 100237 |
| Medit i900/i900c Power Delivery Cable 2.5 m                  | ✓    | ✓     |       | 100238 |
| Medit i900/i900c Power Delivery Cable 2.0 m                  | ✓    | ✓     |       | 100239 |
| Medit i900/i900c/i900M Extended Warranty + 1 year            | ✓    | ✓     | ✓     | 992296 |
| Medit i900/i900c/i900M Large Tip Set (4 pcs)                 | ✓    | ✓     | ✓     | 100240 |
| Medit i900/i900c/i900M Medium Tip Set (4 pcs)                | ✓    | ✓     | ✓     | 100241 |
| Medit i900/i900c/i900M Small Tip Set (4 pcs)                 | ✓    | ✓     | ✓     | 100242 |
| Medit Cart                                                   | ✓    | ✓     | ✓     | 100243 |
| Medit i900M Stand Mount                                      |      |       | ✓     | 100245 |
| Medit i900M Wall Mount Holder                                |      |       | ✓     | 100246 |
| Medit i900M Battery Charger                                  |      |       | ✓     | 100247 |
| Medit i900M Rechargeable Battery (3 pcs.)                    |      |       | ✓     | 100248 |
| Medit i900M Power Adapter                                    |      |       | ✓     | 100249 |
| Medit i900M Charging Cable                                   |      |       | ✓     | 100250 |
| Medit i900M Charging Cable (2m + 30cm)                       |      |       | ✓     | 100250 |

Accessories Medit i700/i700w

| Article                                   | i700w | i700 | REF    |
|-------------------------------------------|-------|------|--------|
| Medit i700/i700w Handpiece Cover          | ✓     | ✓    | 100203 |
| Medit i700w Rechargeable Battery (3 pcs)  | ✓     |      | 100204 |
| Medit i700w Extension Cable Battery       | ✓     |      | 100205 |
| Medit i700w Battery Charger               | ✓     |      | 100206 |
| Medit i700/i700w Normal Tip Set (4 pcs)   | ✓     | ✓    | 100207 |
| Medit i700/i700w Small Tip Set (4 pcs)    | ✓     | ✓    | 100208 |
| Medit i700w Handle                        | ✓     |      | 100209 |
| Medit i700/i700w Calibration Tool         | ✓     | ✓    | 100210 |
| Medit Practice Model                      | ✓     | ✓    | 100211 |
| Medit i700/i700w Wrist Strap              | ✓     | ✓    | 100212 |
| Medit i700/i700w Desktop Cradle           | ✓     | ✓    | 100213 |
| Medit i700w Wall Mount Holder             | ✓     |      | 100214 |
| Medit i700/i700w USB 3.0 Cable            | ✓     | ✓    | 100215 |
| Medit i700w Power Delivery Cable (C to C) | ✓     |      | 100216 |
| Medit i700w Medical Adapter               | ✓     |      | 100217 |
| Medit i700/i700w Power Cord               | ✓     | ✓    | 100218 |
| Medit i700/i700w Monitor Mount            | ✓     | ✓    | 100219 |
| Medit i700w Auto Calibration Tool         | ✓     |      | 100220 |
| Medit i700w Handle-Fit Cradle             | ✓     |      | 100221 |
| Medit i700w 3 m PD Cable                  | ✓     |      | 100222 |
| Medit i700w Multi Hub PKG                 | ✓     |      | 100223 |
| Medit i700w Extended Warranty + 2 years   | ✓     |      | 992281 |
| Medit i700w Extended Warranty + 4 years   | ✓     |      | 992282 |

Accessories Medit i700/i700w

| Article                                | i700w | i700 | REF    |
|----------------------------------------|-------|------|--------|
| Medit i700 Power Hub                   |       | ✓    | 100226 |
| Medit i700 Power Delivery Cable 2m     |       | ✓    | 100227 |
| Medit i700 Power Delivery Cable 3m     |       | ✓    | 100228 |
| Medit i700 Wall Mount Holder           |       | ✓    | 100231 |
| Medit i700 Extended Warranty + 2 years |       | ✓    | 992283 |
| Medit i700 Extended Warranty + 4 years |       | ✓    | 992284 |
| Medit Cart                             | ✓     | ✓    | 100243 |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.



2

**Medit\* Lab Scanner  
& exocad\***

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group



# Lab Scanner and CAD Software

The right choice tailored to your needs

BEGO offers a number of systems combining various scanner models with leading exocad\* CAD software packages for the dental laboratory – tailored to your needs. This allows digitalization of models and conventional impressions as well as design of the corresponding restorations. The scanner and software are assembled and set up for you by our trained technicians. The purchase price includes a training of your technicians in addition to other valuable services.

At BEGO, you'll receive a comprehensive service package with your scanner

|                                                                                  |
|----------------------------------------------------------------------------------|
| Setting up, commissioning, and training                                          |
| Qualified hotline support                                                        |
| Loan unit in the event of damage within 48 h (within the 2-year warranty period) |
| BEGO repair service**                                                            |

You can start working on your first design immediately after installation of the workstation and take advantage of the material and product diversity as well as the outstanding cost effectiveness of BEGO CAD/CAM technology right from the beginning. The coordinated system consisting of a scanner together with scanning and modeling software guarantees you the highest level of quality.

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group  
\*\* Fee-based after expiry of the guarantee  
Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.



# Medit\* T310, T510 & T710

Experience super-fast scanning and high-resolution cameras with high-precision Medit scanners – combined with flexible BEGO exocad software!

Comparison Medit Desktop Scanner

|                      | T710            | T510            | T310            |
|----------------------|-----------------|-----------------|-----------------|
| Auto-elevation       | •               | •               | •               |
| Color texture        | •               | •               | •               |
| Flexible scanning    | •               | •               | optional        |
| Articulator scanning | •               | •               | optional        |
| Replica denture      | •               | optional        | optional        |
| Orthodontic scanning | •               | optional        | optional        |
| Impression scanning  | •               | ×               | ×               |
| Camera Resolution    | mono 5.0 MP × 4 | mono 5.0 MP × 2 | mono 5.0 MP × 2 |
| Accuracy (ISO 12836) | 4 µm            | 7 µm            | 9 µm            |
| Full arch scanspeed  | 8 sec           | 12 sec          | 18 sec          |

Product details

| Scope of delivery                                                                                                                                                                                        | Price / unit |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| Calibration set                                                                                                                                                                                          | inclusive    |
| Standard Module / exocad* DentalCAD <ul style="list-style-type: none"><li>• CAD software for the design of anatomic crowns, copings, bridge frameworks, as well as inlays, onlays, and veneers</li></ul> | inclusive    |
| Virtual Articulator Module <ul style="list-style-type: none"><li>• CAD software allowing users to consider dynamic occlusion</li></ul>                                                                   | inclusive    |


\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.

Medit Scanner incl. software packages

| Medit T310          | REF    |
|---------------------|--------|
| Medit T310 Core     | 280703 |
| Medit T310 Advanced | 280709 |
| Medit T310 Implant  | 280715 |
| Medit T310 Ultimate | 280721 |
|                     |        |
| Medit T510          |        |
| Medit T510 Core     | 280704 |
| Medit T510 Advanced | 280710 |
| Medit T510 Implant  | 280716 |
| Medit T510 Ultimate | 280722 |
|                     |        |
| Medit T710          |        |
| Medit T710 Core     | 280705 |
| Medit T710 Advanced | 280711 |
| Medit T710 Implant  | 280717 |
| Medit T710 Ultimate | 280723 |



\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.



exocad\* software packages

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | REF    |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|
| <b>Software Package Core</b><br>This software package contains the following modules: <ul style="list-style-type: none"><li>DentalCAD Core Version</li><li>Virtual Articulator</li></ul> License update (optional) for 12 months                                                                                                                                                                                                                                                                                                         | 290006 |
| <b>Software Package Advanced</b><br>This software package contains the following modules: <ul style="list-style-type: none"><li>DentalCAD Core Version</li><li>Virtual Articulator</li><li>Provisional Module</li><li>TruSmile</li><li>Tooth Library ZRS Wiedmann</li></ul> License update (optional) for 12 months                                                                                                                                                                                                                      | 992099 |
| <b>Software Package Implant</b><br>This software package contains the following modules: <ul style="list-style-type: none"><li>DentalCAD Core Version</li><li>Virtual Articulator</li><li>Provisional Module</li><li>TruSmile</li><li>Tooth Library ZRS Wiedmann</li><li>Implant Module</li><li>Bar Module</li><li>DICOM Viewer</li></ul> License update (optional) for 12 months                                                                                                                                                        | 290007 |
| <b>Software Package Ultimate</b><br>This software package contains the following modules: <ul style="list-style-type: none"><li>DentalCAD Core Version</li><li>Virtual Articulator</li><li>Provisional Module</li><li>TruSmile</li><li>Tooth Library ZRS Wiedmann</li><li>Implant Module</li><li>Bar Module</li><li>DICOM Viewer</li><li>Model Creator</li><li>Smile Creator</li><li>FullDenture Module</li><li>PartialCAD Module</li><li>Bite Splint Module</li><li>Jaw Motion Import</li></ul> License update (optional) for 12 months | 992100 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 290008 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 992101 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 290009 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 992102 |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

Optional exocad\* software modules

| Implant Module                                                                                                                                                                                                                                                                                                        | REF              |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
| CAD software for the design of abutments and screw-retained bridges<br>License update (optional) for 12 months                                                                                                                                                                                                        | 991979<br>992083 |
| Bar Module                                                                                                                                                                                                                                                                                                            |                  |
| CAD software for the design of standard dental bars as well as advanced custom bar design<br>License update (optional) for 12 months                                                                                                                                                                                  | 991981<br>992091 |
| Model Creator Module                                                                                                                                                                                                                                                                                                  |                  |
| CAD software for the design of models from intra-oral scan data or impression scans<br>License update (optional) for 12 months                                                                                                                                                                                        | 991983<br>992077 |
| Provisional Module                                                                                                                                                                                                                                                                                                    |                  |
| CAD software for the design of temporary crowns and bridges<br>License update (optional) for 12 months                                                                                                                                                                                                                | 991985<br>992073 |
| PartialCAD Module**                                                                                                                                                                                                                                                                                                   |                  |
| CAD software for designing high-quality removable partial denture frameworks<br>License update (optional) for 12 months                                                                                                                                                                                               | 992059<br>992075 |
| Bite Splint Module                                                                                                                                                                                                                                                                                                    |                  |
| CAD software for the design of night guards and individual impression tray<br>License update (optional) for 12 months                                                                                                                                                                                                 | 991987<br>992089 |
| DICOM Viewer Module                                                                                                                                                                                                                                                                                                   |                  |
| Software for visualization of voxel data from CT machines during the design process<br>License update (optional) for 12 months                                                                                                                                                                                        | 991990<br>992087 |
| TruSmile Module                                                                                                                                                                                                                                                                                                       |                  |
| Software for realistic rendering of dental restorations in real-time<br>License update (optional) for 12 months                                                                                                                                                                                                       | 991992<br>992067 |
| Tooth Library Module                                                                                                                                                                                                                                                                                                  |                  |
| Additional tooth library                                                                                                                                                                                                                                                                                              | 991994           |
| Jaw Motion Module                                                                                                                                                                                                                                                                                                     |                  |
| For Import of skull-related survey datan<br>License update (optional) for 12 months                                                                                                                                                                                                                                   | 992026<br>992079 |
| Full Denture Module                                                                                                                                                                                                                                                                                                   |                  |
| CAD software for the design of full denture<br>License update (optional) for 12 months                                                                                                                                                                                                                                | 992028<br>992085 |
| In-CAD Nesting Module                                                                                                                                                                                                                                                                                                 |                  |
| Software supports multi-layered blank and realistic rendering<br>License update (optional) for 12 months                                                                                                                                                                                                              | 992065<br>992081 |
| Smile Design                                                                                                                                                                                                                                                                                                          |                  |
| Software for aesthetic planning using patient photos with contour lines and 3D situations<br>License update (optional) for 12 months                                                                                                                                                                                  | 992063<br>992071 |
| Refresh Fee                                                                                                                                                                                                                                                                                                           |                  |
| Customers without active update contract for their perpetual license now have an opportunity to refresh to the latest software version including purchased modules for a one-time refresh fee plus their annual update contract amount. Program participants will also benefit from additional updates for 12 months. | 992103           |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group  
\*\* Please note that your PC must have at least 4 GB of graphics memory to run PartialCAD!

Optional accessories

|                                       | REF    |
|---------------------------------------|--------|
| Medit* Module Flexible Scanning       | 992127 |
| Medit Module Articulator Scanning     | 100152 |
| Medit Module Replica denture          | 992128 |
| Medit Module Orthodontic scanning     | 992129 |
| Medit Warranty Extension 2 years      | 992130 |
| Medit Calibration panel               | 100153 |
| Medit LED protector                   | 100154 |
| Medit Flexible multi-die jig          | 100155 |
| Medit 16 Single die                   | 100156 |
| Medit Spray supporter                 | 100157 |
| Medit Screw jig                       | 100158 |
| Medit Power cable & adaptor (medical) | 100159 |
| Medit USB 3.0 cable                   | 100160 |
| Medit KAS jig                         | 100161 |
| Medit AM jig                          | 100162 |
| Medit Articulator plate               | 100163 |
| Medit Impression jig                  | 100164 |



\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group  
Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.



# 3Shape – Lab Scanner and CAD Software

You have the choice

Your choice of a new and personal 3Shape system – more flexible than ever before!

- 1

Choose your scanner
- 2

Choose your suitable software package
- 3

Choose the best system for you out of three powerful PCs
- Optional:

Decide on further additional modules and options.
- Always inclusive:

The proven BEGO support and system set-up!

BEGO offers a number of systems combining various scanner models with leading CAD software packages for the dental laboratory – tailored to your needs. This allows digitalization of models and conventional impressions and design of the corresponding restorations. The scanner, PC, and software are assembled and set up for you by our trained technicians. The purchase price includes a training of your technicians in addition to other valuable services.

You can start working on your first design immediately after installation of the workstation and take advantage of the material and product diversity as well as the outstanding cost effectiveness of BEGO CAD/CAM technology right from the beginning. The coordinated system consisting of a scanner as well as scanning and modeling software guarantees you the highest level of quality.

With us you get all scanners with a comprehensive service package

|                                                                                  |   |
|----------------------------------------------------------------------------------|---|
| Setting up, commissioning, and training                                          | ✓ |
| Qualified hotline support                                                        | ✓ |
| Loan unit in the event of damage within 48 h (within the 2-year warranty period) | ✓ |
| BEGO repair service**                                                            | ✓ |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
\*\* Fee-based after expiry of the guarantee  
Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.



# Choose Your 3Shape Scanner – E2, E3, E4, F8

## Step 1

- 2 high-resolution 5.0 MP cameras (E2, E3)
  - 4 high-resolution 5.0 MP cameras (E4, F8)
  - Dual model scanning (F8)
- Integrated multi-die scan technology possible per Add-on
  - Multi line Blue LED
  - Delivery incl. calibration set

Comparison of the 3Shape desktop scanners

|                                            | E2                          | E3                          | E4                          | F8                          |
|--------------------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Number of cameras                          | 2 cameras<br>5.0 megapixels | 2 cameras<br>5.0 megapixels | 4 cameras<br>5.0 megapixels | 2 cameras<br>5.0 megapixels |
| Light technology                           | Multi line<br>Blue LED      | Multi line<br>Blue LED      | Multi line<br>Blue LED      | Multi line<br>Blue LED      |
| Accuracy of crowns and bridges (ISO 12836) | 10 µm                       | 7 µm                        | –                           | –                           |
| Accuracy of implant bars (ISO 12836)       | 12 µm                       | 10 µm                       | 4 µm                        | 4 µm                        |
| Impression scanning                        | ✓ (add-on)                  | ✓ (add-on)                  | ✓                           | ✓                           |
| Multi-Die scanning                         | ✓                           | ✓                           | ✓                           | ✓                           |
| Scanning of textures                       | ✓ (textures)                | ✓ (Color textures)          | ✓ (Color textures)          | ✓ (Color textures)          |
| Dual model scanning                        | –                           | –                           | –                           | ✓                           |
| Full jaw scanning time                     | 24 seconds                  | 18 seconds                  | 9 seconds                   | 9 seconds                   |
| Full jaw impression scanning time          | 72 seconds                  | 64 seconds                  | 45 seconds                  | 45 seconds                  |
| Connection                                 | USB                         | USB                         | USB                         | Ethernet                    |
| REF                                        | 280672                      | 280674                      | 280676                      | 280727                      |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.  
Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.

# Choose Your 3Shape Software Packages

## Step 2

3Shape software packages

|                                                                                                                                                                                                                                                                                                                                                                          | REF       |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| Dental System Scan Only                                                                                                                                                                                                                                                                                                                                                  | inclusive |
| Dental System Crown & Bridge<br>This package includes the following software modules: <ul style="list-style-type: none"><li>• Dental Designer Crown &amp; Bridge</li></ul> Design of full anatomical crowns, bridges, copings, frameworks, Wax-Ups, Impression Scanning, and more                                                                                        | 280684    |
| Dental System Premium<br>This package includes the following software modules: <ul style="list-style-type: none"><li>• Dental Designer Premium</li></ul> Design of full anatomies, bridges, copings, frameworks, and Wax-Ups, splints, digital temps, digital smile design, RealView, telescopes, post & core, Impression Scanning, and more                             | 280682    |
| Dental System Complete Restorative<br>This package includes the following software modules: <ul style="list-style-type: none"><li>• Dental Designer Premium</li><li>• Abutment Designer</li><li>• Implant Bar &amp; Bridge Designer</li><li>• Model Builder All Scans</li><li>• Impression Scanning</li><li>• Removable Partial Design</li><li>• Full Dentures</li></ul> | 280680    |

# Choose PC

## Step 3

3Shape choice PC

|                                                                                                                                                                                 | REF –<br>combined with scanner |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|
| Desktop Standard PC<br>Valid for: <ul style="list-style-type: none"><li>• Scanner E2, Software Crown and Bridge, Premium</li></ul>                                              | 280689                         |
| Desktop High Performance PC<br>Valid for: <ul style="list-style-type: none"><li>• All Scanner, all software packages</li></ul>                                                  | 280690                         |
| Desktop Ultra PC<br>Valid for: <ul style="list-style-type: none"><li>• All Scanner, all software packages</li><li>• Recommended for Implant Bridge and Denture Design</li></ul> | 280691                         |

Note: A monitor is not included.

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

# Choose Your 3Shape LabCare Service

Step 4 (service for your software package)

LabCare service (can be combined with all 3Shape software packages)

| LabCare                                                                          | REF    |
|----------------------------------------------------------------------------------|--------|
| The comprehensive service package that keeps you ahead and constantly productive | 992021 |



\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.

# Choose Your Optional 3Shape Software Modules

Step 5

Optional 3Shape software modules

| 3Shape Abutment Designer                                                                                                                                         | REF              |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
| CAD software for customized bar and bridge constructions<br>Annual license fee as of the 13 <sup>th</sup> month of use                                           | 991916<br>991913 |
| 3Shape Bundle (Abutment Designer and Implant Bar and Bridge Designer)                                                                                            |                  |
| CAD software for customized abutments and CAD software for customized bar and bridge constructions<br>Annual license fee as of the 13 <sup>th</sup> month of use | 991935<br>991936 |
| 3Shape Modellguss Designer                                                                                                                                       |                  |
| CAD software for the Wax-Up of partial denture bases<br>Annual license fee as of the 13 <sup>th</sup> month of use                                               | 991929<br>991930 |
| 3Shape Model Builder – TRIOS                                                                                                                                     |                  |
| CAD-Software for processing intraoral scan data<br>Annual license fee as of the 13 <sup>th</sup> month of use                                                    | 991941<br>-      |
| 3Shape Model Builder ALL Scans                                                                                                                                   |                  |
| Software for processing intraoral scan data<br>Annual license fee as of the 13 <sup>th</sup> month of use                                                        | 991977<br>991921 |
| 3Shape Clear Aligner Studio                                                                                                                                      |                  |
| Complete workflow for orthodontic analysis, diagnosis and treatment planning<br>Annual license fee as of the 13 <sup>th</sup> month of use                       | 992093<br>992095 |

3Shape Full Denture

| CAD software for Full Dentures Modules                                                                                      | REF              |
|-----------------------------------------------------------------------------------------------------------------------------|------------------|
| 3Shape Full Dentures Module <ul style="list-style-type: none"><li>Annual license fee as of the 13th month of use</li></ul>  | 992287<br>992288 |
| CAD software for Full Dentures Modules                                                                                      |                  |
| 3Shape Full Dentures Modules <ul style="list-style-type: none"><li>Annual license fee as of the 13th month of use</li></ul> | 992289<br>992290 |
| CAD software for Full Dentures Modules                                                                                      |                  |
| 3Shape Full Dentures Modules <ul style="list-style-type: none"><li>Annual license fee as of the 13th month of use</li></ul> | 992291<br>992292 |

3Shape software packages

| Dental System Premium                                                                                                                                                                                                                                                                                                                                                                    | REF                                  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|
| You can design an unlimited amount of simple and complex restaurations <ul style="list-style-type: none"><li>Annual license fee as of the 13<sup>th</sup> month of use</li><li>License for another workstation</li><li>Annual license fee as from the second workstation</li></ul>                                                                                                       | 992018<br>991948<br>991949           |
| Dental System Complete Restorative                                                                                                                                                                                                                                                                                                                                                       |                                      |
| You receive unlimited access to every provision case without buying any additional modular units <ul style="list-style-type: none"><li>Upgrade Dental System Premium to Dental System Complete Restorative</li><li>Annual license fee as of the 13<sup>th</sup> month of use</li><li>License for another workstation</li><li>Annual license fee as from the second workstation</li></ul> | 992092<br>992052<br>992123<br>992124 |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

# Choose Optional 3Shape Accessories


## Step 6

Optional 3Shape accessories

| Article                                                                                                                                                                   | REF    |                                                                                       |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|---------------------------------------------------------------------------------------|
| <b>Dental Wax-Up Bridge Fixture</b><br>for fixing complete Wax-Up bridges in the scanner                                                                                  | 100090 |    |
| <b>Scanner Interface Plate</b>                                                                                                                                            | 100092 |    |
| <b>Dental 2 Cast Fixture</b><br>for fixing the occlusal alignment of upper and lower cast during scanning                                                                 | 100093 |    |
| <b>Blu-Tack</b><br>Modelling mass for model fixation                                                                                                                      | 100094 |    |
| <b>Dental Impression Fixture</b><br>for fixing dental impressions in the scanner                                                                                          | 100095 |    |
| <b>Space Explorer 3D-Mouse</b><br>Ergonomically designed 3D mouse for professional users                                                                                  | 100096 |   |
| <b>Transfer Plate for ACR Atex-compatible</b><br>for transfer of the model position from the physical articulator to the dynamic virtual articulator in Dental System     | 100101 |  |
| <b>Transfer Plate for Whip Mix Denar Mark 330</b><br>for transfer of the model position from the physical articulator to the dynamic virtual articulator in Dental System | 100102 |  |
| <b>Transfer Plate for SAM 2P</b><br>for transfer of the model position from the physical articulator to the dynamic virtual articulator in Dental System                  | 100103 |  |
| <b>Transfer Plate for KaVo PROTAR evo</b><br>for transfer of the model position from the physical articulator to the dynamic virtual articulator in Dental System         | 100104 |  |
| <b>Transfer Plate for Generic</b><br>for transfer of the model position from the physical articulator to the dynamic virtual articulator in Dental System                 | 100137 |  |
| <b>Transfer Plate for Shofu MSMP</b><br>for transfer of the model position from the physical articulator to the dynamic virtual articulator in Dental System              | 100142 |  |
| <b>3Shape Artikulator</b><br>for articulation of digital lab models designed in 3Shape's Model Builder                                                                    | 100106 |  |
| <b>Articulator Calibration Object</b><br>for calibrating an articulator                                                                                                   | 100123 |  |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.

Optional 3Shape accessories

| Article                                                                                                                                       | REF    |                                                                                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------|--------|-------------------------------------------------------------------------------------|
| <b>Occlusion Setup Tool incl. 4 Transfer Plates</b><br>allows users to set up upper and lower arches prior to scanning and skip the bite scan | 100124 |  |
| <b>Scan Height Adapter – E scanners</b><br>to adjust the scan height                                                                          | 100133 |  |
| <b>Multi-die Fixture</b><br>for simultaneous scanning of several dies                                                                         | 100134 |  |
| <b>Articulator Holder</b><br>allows the user to mount all 3Shape acknowledge articulators in all E-scanner models                             | 100135 |  |
| <b>3Shape Rubber Bands</b><br>can be used as a replacement or additional Rubber Bands, included in the delivery of Articulator Holder         | 100136 |  |
| <b>CAD Points</b>                                                                                                                             |        |                                                                                     |
| 100 CAD Points                                                                                                                                |        |                                                                                     |
| 200 CAD Points                                                                                                                                |        |                                                                                     |
| 500 CAD Points                                                                                                                                |        |                                                                                     |

3Shape Trade-Up Offers

| Article                                                                           |
|-----------------------------------------------------------------------------------|
| Please consult your local BEGO representative for your individual trade-up offer. |



For a detailed price list for the use of CAD Points, please click here:

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.

4


**Implant and  
Chairside Planning**



# Implant planning

Planning software and BEGO Guide Planning Center

The optimal digital workflow in implantology begins with well thought-out planning. Whether you want to do this yourself and print the template yourself or have our experienced planning team – the decision is entirely up to you. BEGO offers you flexible options to choose from.



exocad\* Implant planning with exoplan

# exocad\* exoplan

| Licences exoplan                                       |  | REF   |
|--------------------------------------------------------|--|-------|
| exocad exoplan License                                 |  | 58621 |
| exocad exoplan License Update                          |  | 99467 |
| exocad exoplan License                                 |  | 58622 |
| exocad exoplan incl. Guide Creator License Update      |  | 99468 |
| exocad exoplan (+2) incl. Guide Creator License        |  | 58566 |
| exocad exoplan (+2) incl. Guide Creator License Update |  | 99469 |
| exocad exoplan (+4) incl. Guide Creator License        |  | 58567 |
| exocad exoplan (+4) incl. Guide Creator License Update |  | 99470 |
| exocad exoplan (+7) incl. Guide Creator License        |  | 58568 |
| exocad exoplan (+7) incl. Guide Creator License Update |  | 99471 |
| exocad exoplan Guide Creator Add-On License            |  | 99465 |
| exocad exoplan Guide Creator Add-On License Update     |  | 99466 |

All exocad exoplan licenses are flex licenses for which an annual license payment is mandatory. The first year is included in the license price.

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.



# 3Shape\* Implant Studio

## Licences Implant Studio

|                                                                      |  | REF   |
|----------------------------------------------------------------------|--|-------|
| 3Shape Implant Studio Stand-Alone                                    |  | 58569 |
| 3Shape Implant Studio Stand-Alone Update                             |  | 99476 |
| 3Shape Implant Studio Stand-Alone Additional Seat                    |  | 58570 |
| 3Shape Implant Studio Stand-Alone Additional Seat Update             |  | 99477 |
| 3Shape Implant Studio Add-On to Dental System                        |  | 99472 |
| 3Shape Implant Studio Add-On to Dental System Update                 |  | 99474 |
| 3Shape Implant Studio Add-On to Dental System Additional Seat        |  | 99473 |
| 3Shape Implant Studio Add-On to Dental System Additional Seat Update |  | 99475 |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

# BEGO Guide Planning Center

Your affordable entry into the world of navigated surgery

Regardless of the implant system used, the BEGO Guide Planning Center enables you to use 3D diagnostics and drilling templates. As a dentist, implantologist and dental technician, you can integrate the system into your daily practice or laboratory routine without having to invest in planning software – with excellent results for you and your patients.

How does it work?

The BEGO Guide Planning Center offers you the opportunity to create a surgical guide, quickly and cost-effectively, without having your own implant planning software. You send the DVT data to our online planning center and we create an implant planning proposal within a very short time.

You will then receive a pilot or fully navigated surgical guide from our digital production facility in Bremen.

Your benefits at a glance

- No investments in software
- Compatibility with all DVT devices
- Delivery of a completely finished drilling template
- Compatible with the BEGO Guide Trays for fully navigated insertion of BEGO Semados® implants
- Pilot surgical guide (sleeve diameter 2.0 mm) for all implant systems

BEGO Guide

|                                                                                                                                                                                                                                                                                                                      |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Planning proposal</b><br>Includes the creation of a planning proposal according to your desired implant positions (regardless of the number of implants to be placed) and the desired implant system;<br>Provision of the planning in the BEGO Guide order portal, incl. online planning discussion, if necessary |
| <b>BEGO Guide (BEGO Semados®)</b><br>Pilot or fully navigated template for BEGO Semados® implants including all sleeves, regardless of the number of implants                                                                                                                                                        |
| <b>BEGO Guide (Foreign implantat)</b><br>Pilot template for third-party implants incl. all sleeves, regardless of the number of implants                                                                                                                                                                             |
| <b>Template design</b><br>Includes the creation of a template file (including scan and matching) either on the basis of the planning from the planning center or by sending a 3Shape* or exoplan planning file (30XZ)                                                                                                |
| <b>Template packages</b><br>Include the planning proposal and the BEGO Guide (BEGO Semados®), as well as the shipping costs (delivery of the template on the next working day by 12.00 noon** ) and the pick up service.                                                                                             |
| <b>Wax-Up</b><br>Design of a digital prosthetic preliminary set-up, price per tooth position                                                                                                                                                                                                                         |
| <b>Pick up service</b><br>Parcel pick up (situation model and CD if applicable) on the desired date, bookable in the BEGO Guide order portal                                                                                                                                                                         |
| <b>Express processing</b><br>For accelerated processing of cases that are released at short notice***                                                                                                                                                                                                                |

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

\*\* After approval of the planning proposal and production of the template (production time usually 2-3 working days)

\*\*\* Express processing shortens the delivery time by an average of one working day

# exocad\* chairsideCAD

You would like to integrate innovative chairside solutions into your practice workflow?

For this purpose we offer exocad's customised modules of the chairsideCAD software.



Licences chairsideCAD

|                                                                                                  | REF            |
|--------------------------------------------------------------------------------------------------|----------------|
| <b>exocad chairsideCAD License</b><br>exocad chairsideCAD Licence Update                         | 58620<br>99451 |
| <b>exocad chairsideCAD Single Visit Bundle</b><br>exocad chairsideCAD Single Visit Bundle Update | 58562<br>99452 |
| <b>exocad chairsideCAD Mill Bundle</b><br>exocad chairsideCAD Mill Bundle Update                 | 58563<br>99453 |
| <b>exocad chairsideCAD Print Bundle</b><br>exocad chairsideCAD Print Bundle Update               | 58564<br>99454 |
| <b>exocad chairsideCAD Ultimate Bundle</b><br>exocad chairsideCAD Ultimate Bundle Update         | 58565<br>99455 |

Modules

|                                                                                                            | REF            |
|------------------------------------------------------------------------------------------------------------|----------------|
| <b>exocad* chairsideCAD Implant Module</b><br>exocad chairsideCAD Implant Module Update                    | 99440<br>99456 |
| <b>exocad chairsideCAD Module Model Creator</b><br>exocad chairsideCAD Module Model Creator Update         | 99441<br>99457 |
| <b>exocad chairsideCAD Bite Splint</b><br>exocad chairsideCAD Bite Splint Update                           | 99442<br>99458 |
| <b>exocad chairsideCAD Provisional Module</b><br>exocad chairsideCAD Provisional Module Update             | 99443<br>99459 |
| <b>exocad chairsideCAD Module Smile Creator</b><br>exocad chairsideCAD Module Smile Creator Update         | 99444<br>99460 |
| <b>exocad chairsideCAD Auto Articulator Module</b><br>exocad chairsideCAD Auto Articulator Module Update   | 99445<br>99461 |
| <b>exocad chairsideCAD Module Jaw Motion Import</b><br>exocad chairsideCAD Module Jaw Motion Import Update | 99446<br>99462 |
| <b>exocad chairsideCAD Module DICOM Viewer</b><br>exocad chairsideCAD Module DICOM Viewer Update           | 99447<br>99463 |
| <b>exocad chairsideCAD TruSmile Module</b><br>exocad chairsideCAD TruSmile Module Update                   | 99448<br>99464 |
| <b>exocad chairsideCAD ZRS Tooth Library</b>                                                               | 99449          |

All chairsideCAD licenses are flex licenses for which an annual license payment is mandatory. The first year is included in the license price.

\* This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

5

Services

5

Services

5

Services

5

Services

# BEGOconnectApp

Real time production status



Follow the real-time status of your BEGO crown, bridge, custom abutment or partial denture orders. The BEGOconnectApp operates as a virtual window into the hightech production center of BEGO and can easily be installed via your BEGO customer support onto your desktop device.


After registering with the BEGO customer number and password, users can check whether their orders have been received by BEGO, are in data preparation or production, have been produced or already sent to the customer. Additionally, the app provides the opportunity to track the delivery status of the courier service provider.

Detailed product information:



# BEGO Education

CAD/CAM par excellence – BEGO CAD/CAM courses for a successful future!



The digitalization of dental processes has resulted in fundamental changes to the dental working environment and demands ever more rapid adaptation to new technologies.

Whether you are a newcomer or an old hand, you can benefit from the know-how of our experienced CAD/CAM specialists and gain an insight into the wide spectrum of new possibilities which this technology offers in terms of materials and production processes.

Within the framework of the various courses you can train in scanning, virtual modeling and construction, amongst others.

We have set up fully equipped workstations and training facilities in Bremen, Germany and at other regional laboratory bases.

Step into the future of customized prosthetics with BEGO Medical!

More information on our entire course program can be found at [www.bego.com](http://www.bego.com).

Detailed product information:





[www.bego.com](http://www.bego.com)

**BEGO Medical GmbH**

Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany  
Tel. +49 421 2028-200 · Fax +49 421 2028-174  
E-Mail [cadcam@bego.com](mailto:cadcam@bego.com) · [www.bego.com](http://www.bego.com)

**BEGO Scan- und Designcenter**

Wilhelm-Herbst-Str. 1 · 28359 Bremen · Germany  
Tel. +49 421 2028-210  
[design@bego.com](mailto:design@bego.com) · [www.bego.com](http://www.bego.com)

Always want to have a look at all BEGO news?

Sign up for our newsletter here: [www.bego.com/newsletter](http://www.bego.com/newsletter)



Errors and changes excepted. Our deliveries of products and services are governed by our Standard Terms of Business (STB), which are available at any time at [www.bego.com](http://www.bego.com) and will be sent upon request. Please note that not all our products and services are available throughout the world, depending on authorities' approvals. For more information, please consult your local BEGO representative.