

# BEGO Interface Overview

Valid from: March 2019

## Interface overview\*

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

<sup>1)</sup> Not available in exocad.

<sup>2)</sup> Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

<sup>3)</sup> In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

<sup>4)</sup> Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection				CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level		
			BEGO Titan Grade 5 Wirobond® M1+	IPS e.max CAD	BEGO Zirkon LT	Katana STML	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT Klebebrücke	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

### Range 0.1 BEGO Implant Systems – Semados® S/RI-Line – Sub-Tec (without Platform Switch)

S 3,25	3,25	■								
S/RI 3,75	3,75	■								
S/RI 4,1	4,1	■	√			√		√	√	√
S/RI 4,5	4,5	■								
S/RI 5,5	5,5	■								

### Range 0.2 BEGO Implant Systems – Semados® Mini (only for bars)

2,7; 2,9; 3,1	3,5	■							√	
---------------	-----	---	--	--	--	--	--	--	---	--

### Range 0.4 BEGO Implant Systems – Semados® Multi<sup>Plus</sup>

3,25; 3,75; 4,1; 4,5; 5,5	Abutment niveau									√
---------------------------	-----------------	--	--	--	--	--	--	--	--	---

### Range 0.5 BEGO Implant Systems – Semados® SC/SCX/RS/R SX/RI – PS-Line (with Platform Switch)

RS/R SX 3,0	3,0	■								
SC/SCX 3,25	3,25	■								
SC/SCX/RS/R SX/RI 3,75	3,75	■	√			√		√	√	√
SC/SCX/RS/R SX/RI 4,1	4,1	■								
SC/SCX/RS/R SX/RI 4,5	4,5	■								
SC/SCX/RS/R SX/RI 5,5	5,5	■								

### Range 1.1 compatible with Straumann synOcta/Standard/Standard Plus/Tapered Effect

3,3; 4,1; 4,8	RN 4,8	■ ■	√			√ <sup>2)</sup>			√	√
6,5	WN 6,5	■								

\* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

**Interface overview\***

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

<sup>1)</sup> Not available in exocad.

<sup>2)</sup> Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

<sup>3)</sup> In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

<sup>4)</sup> Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection					CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level	
			BEGO Titan Grade 5 Wirobond® Ml+	IPS e.max CAD	BEGO Zirkon LT	Katana STML	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT Klebebrücke	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

**Range 1.2 compatible with Straumann BoneLevel**

3,3	NC	■	√				√ <sup>2)</sup>		√	√		
4,1; 4,8	RC	■										

**Range 1.3 compatible with Straumann Narrow Neck CrossFit**

3,5	NNC	■	√				√ <sup>1) 2)</sup>					
-----	-----	---	---	--	--	--	--------------------	--	--	--	--	--

**Range 2.1 compatible with DENTSPLY Implants FRIALIT/XiVe**

3,4	3,4	■										
3,8	3,8	■	√				√			√		
4,5	4,5	■										
5,5	5,5	■										

**Range 2.2 compatible with DENTSPLY Implants Ankylos**

3,5; 4,5; 5,5; 7,0	4,1						√ <sup>1) 2)</sup>					√
--------------------	-----	--	--	--	--	--	--------------------	--	--	--	--	---

**Range 3.1 compatible with Nobel Biocare NobelActive**

3,5	NP	■	√				√		√	√		
4,3; 5,0	RP	■										

**Range 3.2 compatible with Nobel Biocare NobelReplace/Replace Select/NobelSpeedy Replace**

3,5	NP	■										
4,0; 4,3	RP	■	√				√			√		
5,0	WP	■										
6,0	6,0	■										

**Range 3.3 compatible with Nobel Biocare Brånemark System/NobelSpeedy Shorty/Groovy**

3,3	NP	■										
3,75; 4,0	RP	■	√				√			√		
5,0	WP	■										

\* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

**Interface overview\***

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

<sup>1)</sup> Not available in exocad.

<sup>2)</sup> Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

<sup>3)</sup> In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

<sup>4)</sup> Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection					CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level	
			BEGO Titan Grade 5 Wirobond® M+	IPS e.max CAD	BEGO Zirkon LT	Katana STML	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT Klebebrücke	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

**Range 3.4 compatible with Nobel Biocare Multi-unit Abutment**

Nobel Active	Abutment niveau									
Nobel Replace										
Brånemark System										
Astra Tech OsseoSpeed										✓
Straumann synOcta										
Straumann BoneLevel										

**Range 4.1 compatible with BIOMET 3i Osseotite Certain/NanoTite Certain**

3,25; 3/4/3; 4/3	3,4	■								
4,0; 4/5/4; 5/4	4,1	■	✓				✓		✓	
5,0; 5/6/5; 6/5	5,0	■								
6,0; 5/6	6,0	■								

**Range 4.2 compatible with BIOMET 3i Osseotite/NanoTite**

3,25	3,4	■								
3,75; 4,0; 3/4	4,1	■	✓				✓ <sup>1)</sup>		✓	
5,0; 4/5	5,0	■								

**Range 5.1 compatible with DENTSPLY Implants ASTRA TECH OsseoSpeed**

3,5; 4,0	3,5/4,0	■	✓				✓		✓	
4,5; 5,0	4,5/5,0	■								

**Range 5.2 compatible with DENTSPLY Implants ASTRA TECH OsseoSpeed EV**

3,0	3,0	■								
3,6	3,6	■								
4,2	4,2	■	✓			✓		✓	✓	✓
4,8	4,8	■								
5,4	5,4	■								

\* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

**Interface overview\***

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

<sup>1)</sup> Not available in exocad.

<sup>2)</sup> Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

<sup>3)</sup> In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

<sup>4)</sup> Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection					CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level	
			BEGO Titan Grade 5 Wirobond® M+	IPS e.max CAD	BEGO Zirkon LT	Katana STML	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT Klebebrücke	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

**Range 6.1 compatible with Zimmer Dental Tapered Screw-Vent**

3,7; 4,1	3,5											
4,7	4,5		√				√		√		√	
6,0	5,7											

**Range 7.1 compatible with Phibo TSH**

S2 3,3	3,3											
S3 3,6; S4 4,2	4,0		√								√	
S5 4,8	5,0											

**Range 8.1 compatible with CAMLOG CAMLOG SCREW-LINE/ROOT-LINE**

3,3	3,3											
3,8	3,8											
4,3	4,3		√				√ <sup>3)</sup>		√		√	√
5,0	5,0											
6,0	6,0											

**Range 8.2 compatible with CAMLOG CONELOG SCREW-LINE**

3,3	3,3											
3,8	3,8											
4,3	4,3						√ <sup>3)</sup>		√			
5,0	5,0											

**Range 8.3 compatible with CAMLOG iSy**

3,8; 4,4; 5,0	4,5/5,2						√ <sup>1) 2)</sup>					
---------------	---------	--	--	--	--	--	--------------------	--	--	--	--	--

**Range 9.1 compatible with Dyna Dental Dyna Helix**

3,2; 3,6; 4,2; 5,0	3,6/4,0										√	
--------------------	---------	--	--	--	--	--	--	--	--	--	---	--

\* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

**Interface overview\***

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

- <sup>1)</sup> Not available in exocad.
- <sup>2)</sup> Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.
- <sup>3)</sup> In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.
- <sup>4)</sup> Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection					CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level	
			BEGO Titan Grade 5 Wirobond® Ml+	IPS e.max CAD	BEGO Zirkon LT	Katana STML	BEGO PMMA Multicolor	CAD/Cast® alloys	Wirobond® C+	BEGO Zirkon LT Klebebrücke	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

**Range 10.1 compatible with BioComp Implant Level/Abutment Level**

3,4; 4,0	4,0									✓	✓
4,6; 5,2	5,2										

**Range 11.1 compatible with Neoss ProActive**

3,5; 4,0; 4,5; 5,0; 5,5	4,0	■	✓		✓		✓		✓		
-------------------------	-----	---	---	--	---	--	---	--	---	--	--

**Range 14.1 compatible with Dentaurum tioLogic**

3,3	3,4	■									
3,7; 4,2	3,7	■ ■	✓								
4,8; 5,5	4,5	■ ■									

**Range 15.1 compatible with Bredent medical Sky**

3,5	3,5	■									
4,0; 4,5; 5,5	4,0	■ ■ ■	✓		✓		✓				

**Range 16.1 compatible with Medentis ICX-Active-Master**

3,3; 3,75; 4,1; 4,8	4,0	■	✓		✓		✓		✓		
---------------------	-----	---	---	--	---	--	---	--	---	--	--

**Range 17.1 compatible with Global D In-Kone**

4,3	4,3										✓ <sup>4)</sup>
-----	-----	--	--	--	--	--	--	--	--	--	-----------------

**Range 17.2 compatible with Global D EVL**

3,3; 4,0; 5,0	3,3/4,0/5,0	■ ■ ■									✓ <sup>4)</sup>
---------------	-------------	-------	--	--	--	--	--	--	--	--	-----------------

\* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

**BEGO Medical GmbH**  
 Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany  
 Hotline +49 421 2028-200  
 Tel. +49 421 2028-0 · Fax +49 421 2028-174  
 Email info@bego-medical.com · www.bego.com

Always want to have a look at all BEGO news?  
 Sign up for our newsletter here: [www.bego.com/newsletter](http://www.bego.com/newsletter)

