

BEGO Interface Overview

Valid from: May 2019

Interface overview*

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

¹⁾ Not available in exocad.

²⁾ Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

³⁾ In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

⁴⁾ Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection	CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level
			BEGO Titan Grade 5 Wirobond® MI+	IPS e.max CAD BEGO Zirkon LT Katana STML BEGO PMMA Multicolor CAD/Cast® alloys Wirobond® C+	BEGO Zirkon LT Katana STML BEGO PMMA Multicolor Wirobond® C+	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

Range 0.1 BEGO Implant Systems – Semados® S/RI-Line – Sub-Tec (without Platform Switch)

S 3,25	3,25	■					
S/RI 3,75	3,75	■					
S/RI 4,1	4,1	■	√		√	√	√
S/RI 4,5	4,5	■					
S/RI 5,5	5,5	■					

Range 0.2 BEGO Implant Systems – Semados® Mini (only for bars)

2,7; 2,9; 3,1	3,5	■				√	
---------------	-----	---	--	--	--	---	--

Range 0.4 BEGO Implant Systems – Semados® Multi^{Plus}

3,25; 3,75; 4,1; 4,5; 5,5	Abutment Niveau						√
---------------------------	-----------------	--	--	--	--	--	---

Range 0.5 BEGO Implant Systems – Semados® SC/SCX/RS/R SX/RI – PS-Line (with Platform Switch)

RS/R SX 3,0	3,0	■					
SC/SCX 3,25	3,25	■					
SC/SCX/RS/R SX/RI 3,75	3,75	■	√		√	√	√
SC/SCX/RS/R SX/RI 4,1	4,1	■				√	√
SC/SCX/RS/R SX/RI 4,5	4,5	■					
SC/SCX/RS/R SX/RI 5,5	5,5	■					

Range 1.1 compatible with Straumann synOcta/Standard/Standard Plus/Tapered Effect (Variobase Abutment)

3,3; 4,1; 4,8	RN 4,8	■ ■	√		√ ²⁾		√
6,5	WN 6,5	■					√

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.

Interface overview*

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

¹⁾ Not available in exocad.

²⁾ Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

³⁾ In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

⁴⁾ Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection	CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level
			BEGO Titan Grade 5 Wirobond® M+	IPS e.max CAD BEGO Zirkon LT Katana STML BEGO PMMA Multicolor CAD/Cast® alloys Wirobond® C+	BEGO Zirkon LT Katana STML BEGO PMMA Multicolor Wirobond® C+	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

Range 1.2 compatible with Straumann BoneLevel

3,3	NC	■	✓		✓ ²⁾	✓	
4,1; 4,8	RC	■					

Range 1.3 compatible with Straumann Narrow Neck CrossFit

3,5	NNC	■	✓		✓ ^{1) 2)}		
-----	-----	---	---	--	--------------------	--	--

Range 2.1 compatible with DENTSPLY Implants FRIALIT/XiVE

3,4	3,4	■					
3,8	3,8	■	✓		✓		
4,5	4,5	■				✓	
5,5	5,5	■					

Range 2.2 compatible with DENTSPLY Implants Ankylos (Ankylos CX Ti-Interface)

3,5; 4,5; 5,5; 7,0	4,1				✓ ^{1) 2)}		✓
--------------------	-----	--	--	--	--------------------	--	---

Range 3.1 compatible with Nobel Biocare NobelActive

3,5	NP	■	✓		✓	✓	
4,3; 5,0	RP	■					

Range 3.2 compatible with Nobel Biocare NobelReplace/Replace Select/NobelSpeedy Replace

3,5	NP	■					
4,0; 4,3	RP	■	✓		✓		
5,0	WP	■				✓	
6,0	6,0	■					

Range 3.3 compatible with Nobel Biocare Brånemark System/NobelSpeedy Shorty/Groovy

3,3	NP	■					
3,75; 4,0	RP	■	✓		✓		✓
5,0	WP	■					

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.

Interface overview*

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

- ¹⁾ Not available in exocad.
- ²⁾ Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.
- ³⁾ In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.
- ⁴⁾ Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection	CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level
			BEGO Titan Grade 5 Wirobond® M+	IPS e.max CAD BEGO Zirkon LT Katana STML BEGO PMMA Multicolor CAD/Cast® alloys Wirobond® C+	BEGO Zirkon LT Katana STML BEGO PMMA Multicolor Wirobond® C+	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+
							

Range 3.4 compatible with Nobel Biocare Multi-unit Abutment

Nobel Active	Abutment Niveau						✓
Nobel Replace							
Brånemark System							
Astra Tech OsseoSpeed							
Straumann synOcta							
Straumann BoneLevel							

Range 4.1 compatible with BIOMET 3i Osseotite Certain/NanoTite Certain

3,25; 3/4/3; 4/3	3,4	■						
4,0; 4/5/4; 5/4	4,1	■	✓		✓		✓	
5,0; 5/6/5; 6/5	5,0	■						
6,0; 5/6	6,0	■						

Range 4.2 compatible with BIOMET 3i Osseotite/NanoTite

3,25	3,4	■					✓	
3,75; 4,0; 3/4	4,1	■	✓		✓ ¹⁾		✓	
5,0; 4/5	5,0	■						

Range 5.1 compatible with DENTSPLY Implants ASTRA TECH OsseoSpeed

3,5; 4,0	3,5/4,0	■					✓	
4,5; 5,0	4,5/5,0	■	✓		✓			

Range 5.2 compatible with DENTSPLY Implants ASTRA TECH OsseoSpeed EV

3,0	3,0	■						
3,6	3,6	■						
4,2	4,2	■	✓		✓		✓	✓
4,8	4,8	■						
5,4	5,4	■						

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.

Interface overview*

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

¹⁾ Not available in exocad.

²⁾ Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

³⁾ In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

⁴⁾ Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection	CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level
			BEGO Titan Grade 5 Wirobond® M+	IPS e.max CAD BEGO Zirkon LT Katana STML BEGO PMMA Multicolor CAD/Cast® alloys Wirobond® C+	BEGO Zirkon LT Katana STML BEGO PMMA Multicolor Wirobond® C+	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

Range 6.1 compatible with Zimmer Dental Tapered Screw-Vent

3,7; 4,1	3,5	■					
4,7	4,5	■	✓		✓	✓	
6,0	5,7	■					

Range 7.1 compatible with Phibo TSH

S2 3,3	3,3	■					
S3 3,6; S4 4,2	4,0	■	✓			✓	
S5 4,8	5,0	■					

Range 8.1 compatible with CAMLOG CAMLOG SCREW-LINE/ROOT-LINE

3,3	3,3	■					
3,8	3,8	■					
4,3	4,3	■	✓		✓ ³⁾	✓	✓
5,0	5,0	■					
6,0	6,0	■					

Range 8.2 compatible with CAMLOG CONELOG SCREW-LINE

3,3	3,3	■					
3,8	3,8	■					
4,3	4,3	■			✓ ³⁾	✓	
5,0	5,0	■					

Range 8.3 compatible with CAMLOG iSy (iSy Hybrid Abutment)

3,8; 4,4; 5,0	4,5/5,2	■			✓ ^{1) 2)}		
---------------	---------	---	--	--	--------------------	--	--

Range 9.1 compatible with Dyna Dental Dyna Helix

3,2; 3,6; 4,2; 5,0	3,6/4,0	■				✓	
--------------------	---------	---	--	--	--	---	--

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.

Interface overview*

Unless specified otherwise, the options listed for the corresponding implant systems are available for exocad and 3Shape users.

¹⁾ Not available in exocad.

²⁾ Two-piece abutment on original adhesive abutments from the manufacturer. Please procure the necessary components from the respective manufacturer.

³⁾ In contrast to the nt-trading titanium base, the original adhesive abutment is supplied.

⁴⁾ Please procure the necessary components from the respective manufacturer.

Implant diameter	Platform	Color coding	CADAbut Full	CADAbut Duo with rotation protection	CADAbut Duo without rotation protection	CADBase implant level	CADBase abutment level
			BEGO Titan Grade 5 Wirobond® M+	IPS e.max CAD BEGO Zirkon LT Katana STML BEGO PMMA Multicolor CAD/Cast® alloys Wirobond® C+	BEGO Zirkon LT Katana STML BEGO PMMA Multicolor Wirobond® C+	BEGO Titan Grade 5 Wirobond® M+	BEGO Titan Grade 5 Wirobond® M+

Range 10.1 compatible with BioComp Implant Level/Abutment Level

3,4; 4,0	4,0						
4,6; 5,2	5,2					✓	✓

Range 11.1 compatible with Neoss ProActive

3,5; 4,0; 4,5; 5,0; 5,5	4,0	■	✓		✓		✓
-------------------------	-----	---	---	--	---	--	---

Range 13.1 compatible with BioHorizons Internal

3,0	3,0	■					
3,5	3,5	■					
4,5	4,5	■		✓			✓
5,7	5,7	■					

Range 14.1 compatible with Dentaurem tioLogic

3,3	3,4	■					
3,7; 4,2	3,7	■	✓				
4,8; 5,5	4,5	■					

Range 15.1 compatible with Bredent medical Sky

3,5	3,5	■	✓				
4,0; 4,5; 5,5	4,0	■		✓			✓

Range 16.1 compatible with Medentis ICX-Active-Master

3,75; 4,1; 4,8	4,0	■	✓		✓		✓
----------------	-----	---	---	--	---	--	---

Range 17.1 compatible with Global D In-Kone

4,3	4,3						✓ ⁴⁾
-----	-----	--	--	--	--	--	-----------------

Serie 17.2 compatible with Global D EVL

3,3; 4,0; 5,0	3,3/4,0/ 5,0	■					✓ ⁴⁾
---------------	--------------	---	--	--	--	--	-----------------

* The following symbols are commercial designations/registered trademarks of companies which are not part of the BEGO company group, with the exception of BEGO Implant Systems (Semados®).

Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.