

PRESS RELEASE

Bremen, 2nd of May, 2018

Osteotome Sets for BEGO Semados[®] implants

The Osteotome Sets for BEGO Semados[®] implants are now provided by Ustomed * and are of high quality and easy to handle.

The osteotomes from the Osteotome Sets are used for the condensing preparation of the local bone in case of soft bone conditions (D3 / D4, according to Lekholm & Zarb) to achieve a sufficient primary stability upon the insertion of the implants. They preserve the bone, that would be removed during an implant bed preparation with rotating tools. The osteotomes are also used to lift the maxillary sinus floor. The depth marks, which are adjusted to the implant lengths, enables a safe preparation of the cavity.

To order the Osteotomes Sets visit the Ustomed* at www.ustomed.de/en under Bego_RS RSX_Set and Bego_SC SCX_Set.

Image caption (©Ustomed) Osteotome-Set for BEGO Semados[®]

RS/RSX Implants <https://share.bego.com/s/LyJ5zyrQsdX57e3>

* This symbol is a commercial name / registered trademark of a business that does not belong to the BEGO Group.

Press contact: *Maike Wachendorf, Product Marketing*

Phone +49 421 2028-267 | Fax: +49 421 2028-44267 |

[*wachendorf@bego.com*](mailto:wachendorf@bego.com)

About BEGO Implant Systems

Bremen-based BEGO Implant Systems GmbH & Co. KG is a thriving, growth-oriented company in the dental implant industry. The company has been developing and manufacturing dental implants and accessories for implant-based treatment of patients around the world since 1990. Dental implants "Made by BEGO" epitomise top German-made products at a fair price, delivering a perfect combination of safety, durability, aesthetics and reliability. BEGO Implant Systems GmbH & Co. KG has patented many of its developments.